GENERAL COMMANDS

!1	Repeat the previous command
!IOFF	Disable timed homing
!IDON	Enable timed homing
INIT or IHOME	Initialize (home) all axes
ISERN	Show the serial number, returns: "-SERN: xxxx"
POSI	Display position of all axes
!STAT	Display temperatures, vacuum, pressure, status of bottles
!REVS	Report the software revisions

PLATE X/Y CONTROL

	Move plate X-axis to home
	Move plate Y-axis to nome
	Move both plate X & Y to home
	Nove plate X-axis to coordinate xxxx
	Move plate Y-axis to coordinate yyyy
	Move plate X & Y to coordinates
	Reference the plate Y-axis again the indexing sensor
	Move plate X, time the plate moving in the X axis for N number of times
IPLYT	Move plate Y, time the plate moving in the Y axis for N number of times

MIXING

MXONtttt	Turn on plate mixer, tttt = time in seconds
MXOF	Turn off mixer

WASH HEAD CONTROL

!WSHH	Move wash head to home
!WSHMzzzz	Move wash head to coordinate zzzz

PROBE X/Z CONTROL

PRXH	Move probe X-axis to home
PRZH	Move probe Z-axis to home
PRZMxxxx	Move probe X to coordinate xxxx
PRZMzzzz	Move probe Z to coordinate zzzz
PRBMxxxxzzz	Move probe X & Z to coordinates.
XTIM	Move probe X, time the probe moving in the X axis for N number of times.
ZTIM	Move probe Z, time the probe moving in the Z axis for N number of times.
PFFR	Find fluid level on rack (sense liquid level)
PFFR	Hind fluid level on rack (sense liquid level)
PARK	Move probe to waste area of wash cup.

RACK CONTROL

!RK1H	Move rack 1 (left) to home
!RK2H	Move rack 2 (right) to home
!RK1Myyyy	Move rack 1 to coordinate yyyy
!RK2Myyyy	Move rack 2 to coordinate yyyy
!RKSMxxxxyyyy	Move both racks at the same time. xxxx = coord.for R1, yyyy = coord. for R2
!R1TM	Move rack 1, time rack 1 moving in the Y axis for N number of times.
!R2TM	Move rack 2, time rack 2 moving in the Y axis for N number or times.

PRESSURE/VACUUM AND VALVE CONTROL

PRON	Pressure system on
VAON	Vacuum system on
VOFF	Vacuum system off
STBY	Standby. Turns off both pressure and vacuum.
!VALT	Check washer mechanism valves, this command will activate them one by one several times
	•

DILUTER CONTROL

IVAL0 or IVAL1	Turn/switch valve to position 0 or 1, use this to determine if the valve is working, a click should be
heard.	
!VALR	Set valve type to solenoid, use this only if valve has lost configuration.
!VALS	Set valve type to rotary, use this only if valve has lost configuration.
!WPRBvvvv	Wash probe tip. vvvv = volume in uL, 2000 = 2.0mL, if vvvv is not specified the default is used.
!DINI	Initialize the syringe pumps
!DPRIn	Prime the syringe pumps. n = number of times to prime.
DASPvvvvSx	Aspirate into syringe. vvvv = volume in uL (large or small chosen automatically). x = syringe speed.
!DDSPvvvvSx	Dispense from the syringe. vvvv = volume in uL, x = syringe speed.

PARAMETERS

PARM	Display the current parameters
!PARMnapppp	Edit a parameter. n = par.#, a = par. label (xyz,ect.), pppp = new par. Value
MSPDssss	Set mix speed duty cycle, ssss = duty (0030 is the default). A larger number results in a slower mix speed.
MAXZxxxx	Max probe depth, xxxx = maximum depth over rack 1
MXZ2xxxx	Max probe depth, xxxx = maximum depth over rack 2

PHOTOMETER

IFLTS	Display the available wavelengths.
ILPON	Turn the lamps on.
ILOFF	Turn the lamps off.
!RFLI !RWELxnnpdnB !VWELxnn !RAIR !BLNKxnn	Bisplay the filter voltages of the 4 channels at all 8 wavelengths. Read the specified well. $x = strip(A-H)$, $nn = strip(1-12)$, $pd = filters$, $n = #$ of readings. Display the filter volts at well position. $x = strip(A-H)$, $nn = strip(1-12)$ Take and store an air reading. Read and store blank values, $x = strip(A-H)$, $nn = strip(1-12)$, entered for channel #2 ((right rear) (all 4 channels, all 8 wavelengths, refer to error 523.))

WASHER FUNCTIONS

!WPRI	Prime the wash system one time.
IRPRI	Prime the rinse system one time.
!ASPSsseed	Aspirate on the plate. ss = start strip(1-12),ee = end strip(1-12),d = "D" for double
IDSPSsseevvv	Dispense on the plate. ss = start strip(1-12),ee = end strip(1-12),vvv = vol. in uL.
!ASDSsseevvv	Do an aspirate/dispense on plate. ss = start strip(1-12), ee = end strip(1-12), vvv = vol. in uL.
!RWSHnnvvvx	Wash row nn, vvv vol in uL, x = number of cycles.
ICASP	Check wash head aspiration (must change 0.5" open versus closed system)

PROBE MACROS

!GOWLxnn	Go to the specified well. $x = strip(A-H)$, nn = strip(1-12).
GOWDxnnvvvvSs	Go to the specified well and dispense, x=strip(A-H), nn =strip(1-12), vvvv = vol.in uL, s = speed(0-5)
!GOR1nn	Go to the specified position in Rack 1, nn = position number.
!GOR2nn	Go to the specified position in Rack 2, nn = position number.
!GR1AnnvvvvSs	Go to the position in Rack 1 and aspirate, nn = position, vvvv = vol. in uL, s = speed (0-5)
!GR1DnnvvvvSs	Go to the position in Rack 1 and dispense, nn = position, vvvv = vol. in uL, s = speed (0-5)
!GR2AnnvvvvSs	Go to the position in Rack 2 and aspirate, nn = position, vvvv = vol. in uL, s = speed (0-5)
!GR2DnnvvvvSs	Go to the position in Rack 2 and dispense, nn = position, vvvv = vol. in uL, s = speed (0-5)