

### 7.3.2 Calibration Errors

<b>System Leak Test</b>		
Error 6 – Pressure Drop	Excessive Pressure Drop During System Test	<ol style="list-style-type: none"> <li>1. Ensure closed patient system.</li> <li>2. Run system test with single tubing limb connect between outlet &amp; exhalation cover.</li> <li>3. Remove system tubing and verify flow from PV1</li> <li>4. Perform internal leakage test.</li> <li>5. Replace Sensor PCB</li> </ol>
Error 7 – Pressure Rise	Excessive Pressure Increase During System Test	<ol style="list-style-type: none"> <li>1. Forward leak through PV1 valve.</li> <li>2. Forward leak through Solenoid 4 if nebulizer is connected.</li> <li>3. Replace Sensor PCB.</li> </ol>
Error 8 – Max Time To Pressure	Unable To Adequately Pressurize Patient System	<ol style="list-style-type: none"> <li>1. Ensure closed patient system.</li> <li>2. Run system test with single tubing limb connect between outlet &amp; exhalation cover.</li> <li>3. Remove system tubing and verify flow from PV1</li> <li>4. Perform internal leakage test.</li> <li>5. Replace Sensor PCB</li> </ol>
Error 10 – Deviation High	Flow At Transducer dP2 Out Of Range	<ol style="list-style-type: none"> <li>1. Verify internal flow sensor calibration, recalibrate as necessary.</li> </ol>
Error 11 – Error Emptying Tank	Tank Pressure Remains Above 100mbar	<ol style="list-style-type: none"> <li>1. Ensure that patient wye is open at the beginning of system test.</li> <li>2. Verify that PV1 is opening and flow is evident from to patient port.</li> <li>3. Replace Inspiratory Valve PV1</li> <li>4. Verify no leakage through blender system.</li> </ol>
<b>Flow Sensor Calibration</b>		
Error 12 – High Pressure	Pressure Out Of Range	<ol style="list-style-type: none"> <li>1. Ensure that the proximal sensor is connected to the tubing system and that its outlet is open.</li> <li>2. Replace proximal sensor</li> <li>3. Verify internal flow sensor calibration, recalibrate as necessary.</li> <li>4. Verify PV1 calibration, recalibrate as necessary.</li> </ol>
Error 13 – Deviation High	Zero Adjustment dP2	<ol style="list-style-type: none"> <li>1. Perform re-zero of pressure transducers, fab test 5</li> <li>2. Replace Sensor PCB.</li> </ol>
Error 14 – Deviation High	Zero Adjustment dP3	<ol style="list-style-type: none"> <li>1. Perform re-zero of pressure transducers, fab test 5</li> <li>2. Replace Sensor PCB.</li> </ol>
Error 15 – Deviation High	Zero Adjustment P3	<ol style="list-style-type: none"> <li>1. Perform re-zero of pressure transducers, fab test 5</li> <li>2. Replace Sensor PCB.</li> </ol>

Error 16 – Deviation High	Differential pressure too low @ 30 lpm	<ol style="list-style-type: none"> <li>1. Ensure that system is leak free by running system test.</li> <li>2. Ensure that the proximal sensor is correctly connected to the tubing system and that its outlet is open.</li> <li>3. Replace proximal sensor</li> <li>4. Perform re-zero of pressure transducers, fab test 5</li> <li>5. Verify internal flow sensor calibration, recalibrate as necessary.</li> <li>6. Replace Sensor PCB</li> </ol>
Error 17 – Deviation High	Differential pressure between Adult / Infant limits @ 30 lpm	<ol style="list-style-type: none"> <li>1. Ensure that system is leak free by running system test.</li> <li>2. Ensure that the proximal sensor is correctly connected to the tubing system and that its outlet is open.</li> <li>3. Replace proximal sensor</li> <li>4. Perform re-zero of pressure transducers, fab test 5</li> <li>5. Verify internal flow sensor calibration, recalibrate as necessary.</li> <li>6. Replace Sensor PCB</li> </ol>
Error 18 – Deviation High	Differential pressure too high @ 30 lpm	<ol style="list-style-type: none"> <li>1. Ensure that system is leak free by running system test.</li> <li>2. Ensure that the proximal sensor is correctly connected to the tubing system and that its outlet is open.</li> <li>3. Replace proximal sensor</li> <li>4. Perform re-zero of pressure transducers, fab test 5</li> <li>5. Verify internal flow sensor calibration, recalibrate as necessary.</li> <li>6. Replace Sensor PCB</li> </ol>
Error 19 – Deviation High	Differential pressure too low @ low flow level	<ol style="list-style-type: none"> <li>1. Ensure that system is leak free by running system test.</li> <li>2. Ensure that the proximal sensor is correctly connected to the tubing system and that its outlet is open.</li> <li>3. Replace proximal sensor</li> <li>4. Perform re-zero of pressure transducers, fab test 5</li> <li>5. Verify internal flow sensor calibration, recalibrate as necessary.</li> <li>6. Replace Sensor PCB</li> </ol>

Error 20 – Deviation High	Differential pressure too high @ low flow level	<ol style="list-style-type: none"> <li>1. Ensure that system is leak free by running system test.</li> <li>2. Ensure that the proximal sensor is correctly connected to the tubing system and that its outlet is open.</li> <li>3. Replace proximal sensor</li> <li>4. Perform re-zero of pressure transducers, fab test 5</li> <li>5. Verify internal flow sensor calibration, recalibrate as necessary.</li> <li>6. Replace Sensor PCB</li> </ol>
Error 21 – Deviation High	Coefficient (a) too low	<ol style="list-style-type: none"> <li>1. Ensure that system is leak free by running system test.</li> <li>2. Ensure that the proximal sensor is correctly connected to the tubing system and that its outlet is open.</li> <li>3. Replace proximal sensor</li> <li>4. Perform re-zero of pressure transducers, fab test 5</li> <li>5. Verify internal flow sensor calibration, recalibrate as necessary.</li> <li>6. Replace Sensor PCB</li> </ol>
Error 22 – Deviation High	Coefficient (a) too high	<ol style="list-style-type: none"> <li>1. Ensure that system is leak free by running system test.</li> <li>2. Ensure that the proximal sensor is correctly connected to the tubing system and that its outlet is open.</li> <li>3. Replace proximal sensor</li> <li>4. Perform re-zero of pressure transducers, fab test 5</li> <li>5. Verify internal flow sensor calibration, recalibrate as necessary.</li> <li>6. Replace Sensor PCB</li> </ol>
Error 23 – Deviation High	Coefficient (b) too low	<ol style="list-style-type: none"> <li>1. Ensure that system is leak free by running system test.</li> <li>2. Ensure that the proximal sensor is correctly connected to the tubing system and that its outlet is open.</li> <li>3. Replace proximal sensor</li> <li>4. Perform re-zero of pressure transducers, fab test 5</li> <li>5. Verify internal flow sensor calibration, recalibrate as necessary.</li> <li>6. Replace Sensor PCB</li> </ol>

Error 24 – Deviation High	Coefficient (b) too high	<ol style="list-style-type: none"> <li>1. Ensure that system is leak free by running system test.</li> <li>2. Ensure that the proximal sensor is correctly connected to the tubing system and that its outlet is open.</li> <li>3. Replace proximal sensor</li> <li>4. Perform re-zero of pressure transducers, fab test 5</li> <li>5. Verify internal flow sensor calibration, recalibrate as necessary.</li> <li>6. Replace Sensor PCB</li> </ol>
Error 25 – Saving Data	NVRAM Damaged	<ol style="list-style-type: none"> <li>1 Replace proximal sensor</li> <li>2.Perform NVRAM test</li> <li>3.Clear &amp; Test NVRAM, recalibrate device.</li> <li>4. Replace controller PCB</li> </ol>
Error 26 – Emptying Tank	Error Emptying Tank	<ol style="list-style-type: none"> <li>1. Ensure that patient wye is open at the beginning of system test.</li> <li>2. Verify that PV1 is opening and flow is evident from to patient port.</li> <li>3. Replace Inspiratory Valve PV1</li> <li>4. Verify no leakage through blender system.</li> </ol>
Error 27 – Low Pressure	Pressure P2 < 15 mbar	<ol style="list-style-type: none"> <li>1. Ensure that that the flow sensor outlet is blocked during compensation test.</li> <li>2. Replace proximal flow sensor</li> <li>3. Perform calibration of proportion valve PV1.</li> </ol>
Error 28 – High Pressure	Pressure P2 > 40 mbar	<ol style="list-style-type: none"> <li>1. Replace proximal flow sensor</li> <li>2. Perform calibration of proportion valve PV1.</li> </ol>
<b>Oxygen Sensor</b>		
Error 2 – Saving Data	NVRAM Damaged.	<ol style="list-style-type: none"> <li>1.Replace O2 Sensor</li> <li>2.Perform NVRAM test</li> <li>3.Clear &amp; Test NVRAM, recalibrate device.</li> <li>4. Replace controller PCB</li> </ol>
Error 3 – Deviation High	Error during calibration at 100% setting.	<ol style="list-style-type: none"> <li>1. Ensure O2 supply is connected and adequate flow available.</li> <li>2. Ensure O2 Sensor connected correctly.</li> <li>3. Replace O2 Sensor.</li> <li>4. Confirm that inlet check valves CV1/2 are correctly installed and undamaged, replace as necessary.</li> <li>5. Replace Sensor PCB.</li> <li>6. Replace O2 Sensor interface block.</li> </ol>