STORZ

KARL STORZ — ENDOSKOPE

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3. Controls, displays, connectors, and their uses

- ① Power switch ("O" = off)
- ② Mode of operation indicator lamps
- ③ Flow-rate indicator
- ④ Flow-rate setpoint indicator
- ⑤ Irrigation pressure level indicator
- ⑥ Irrigation pressure limit indicator
- ⑦ Suction pressure limit indicator
- ⑧ Suction pressure setpoint indicator
- ⑨ Pushbutton, releases pump rollers
- ⑩ Label strip (replaceable)
- ⑪ Pressure transducer
- ⑫ Pilot lamp
- ⑬ Irrigation START/STOP button
- ⑭ Suction line connector
- ⑮ Suction START/STOP button
- ⑯ ± buttons for adjusting suction pressure setpoint
- ⑰ ± buttons for adjusting irrigation pressure setpoint
- ⑱ ± buttons for adjusting flow-rate setpoint
- ⑲ Footswitch connector
- ⑳ SCB* connectors
- ㉑ Potential equalization connector
- ㉒ Power fuses
- ㉓ Power supply cord receptacle
- ㉔ Fluid collection bottle
- ㉕ Fluid collection bottle vacuum line fitting (with bacterial filter)
- ㉖ Fluid collection bottle instrument line fitting

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
* The KARL STORZ SCB interface (KARL STORZ Communication Bus), based on the CAN field bus, permits remote control of equipment functions, as well as remote display of equipment parameters.

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
4. Symbols employed

	Follow instructions for use
	Read the instructions carefully before operating the equipment!
	Flow
	Irrigation pressure
	Suction pressure
	Suction inlet fitting
	Pressure transducer
	START/STOP
	Potential equalization connector
	Footswitch connector
	Applied part of type BF
	Alternating current
	Electronic information product pollution control (China RoHS)
	This device has been labelled in Accordance with the European Directive on waste electrical and electronic equipment (WEEE).
	Manufacturer

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1. Important information for users of KARL STORZ devices	III
2. Images of the unit	V
3. Controls, displays, connectors, and their uses	VII
4. Symbols employed	VIII
5. General information	3
5.1 Device description	3
5.2 Property rights	3
6. Safety instructions	4
6.1 Explanation of warnings and cautions	4
6.2 Intended use	8
6.3 User qualification	9
6.4 Training in the operation and function of the device	9
6.5 Patient profile	9
6.6 Intended region of treatment on the body	9
6.7 User profile of physician and assistants	10
6.8 Intended conditions of use	10
6.9 Safety precautions at the site of installation	11
6.10 Safety precautions when operating the unit	11
6.11 Safety features	12
7. Installation and operating instructions	13
7.1 Unpacking the equipment	13
7.2 Basic equipment	13
7.3 Installation	13
7.4 KARL STORZ SCB	15
7.5 Mode of operation selection	15
7.6 Operating the unit	16
7.7 Set value confirmation	16
7.8 Power-up with irrigation START/STOP button pressed down	16
7.9 Inserting the tubes	17
7.10 Removing the tubing set	18
7.11 Connecting the tubing set for suction (optional)	19
7.12 Presetting irrigation parameters	20

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7.13	Presetting suction pressure.....	20
7.14	Test for proper operation and air-relief of the tubing system.....	21
7.15	Initiating irrigation	22
7.16	Initiating suction	22
8.	Maintenance	23
8.1	Fuse replacement.....	23
8.2	Reprocessing	24
8.3	Maintenance and safety check	31
8.4	Servicing and repair	32
8.5	Disposal	32
8.6	Repair program	33
8.7	Important information	33
8.8	Limitation of liability	34
8.9	Manufacturer's warranty.....	34
9.	Technical description	35
9.1	Troubleshooting	35
9.2	Error codes	37
9.3	Technical data.....	38
9.4	Technical documentation	39
9.5	General circuit diagram	40
10.	Spare parts, recommended accessories	41
10.1	List of Spare parts.....	41
10.2	Recommended accessories	41
11.	Electromagnetic Compatibility (EMC) Information	42
12.	Subsidiaries	58

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5. General information

5.1 Device description

The HAMOU® ENDOMAT® is a combination irrigation/suction pump for use in gynecological and laparoscopic surgery. Use of the most advanced technologies allows its application in both hysteroscopic and laparoscopic surgery. Adaptation to the type of surgery intended via provision of the relevant optimal operating parameters occurs automatically when the correct type of tubing set for the hysteroscopic or laparoscopic application is inserted. Inserting pressure lines has been greatly simplified through use of a redesigned reel pump. The unit electronically monitors the correct positioning of pressure lines.

Operator adjustments and unit operation are made using precision pushbutton switches. Displays for setpoints and current actual performance data, which are positioned next to each other, allow optimal monitoring of unit operational status. Electronic safety circuitry interrupts irrigation/suction operation in the event that unit performance persistently deviates from the setpoint, and issues an audible alarm to warn the operator. An electronic self-check routine that checks all unit components each time the unit is switched on and notifies the operator of any failures detected further enhances the HAMOU® ENDOMAT®'s high standards of safe operation.

5.2 Property rights

This product is protected in the USA by (at least one of the following) US Patent No(s) 5,788,688; 6,397,286; 6,484,221; 6,824,539.

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6. Safety instructions

Safety instructions are measures intended to protect the user and patients from the risks which could arise through the use of the system.

6.1 Explanation of warnings and cautions

Please read this manual and follow its instructions carefully. The words **Warning**, **Caution**, and **Note** convey special meanings. Wherever they are used in this manual, they should be carefully reviewed to ensure the safe and effective operation of this product. To make the words stand out more clearly, they are accompanied by a pictogram.

WARNING: A Warning indicates that the personal safety of the patient or physician may be involved. Disregarding a Warning could result in injury to the patient or physician.

CAUTION: A Caution indicates that particular service procedures or precautions must be followed to avoid possible damage to the product.

NOTE: A Note indicates special information about operating the product, or clarifies important information.

WARNING: Read this instruction manual thoroughly and be familiar with its contents prior to using this equipment. Before using the unit, read the following safety instructions carefully to avoid putting your patients, personnel or yourself at risk.

NOTE: Any damage to the unit resulting from incorrect operation is not covered by the manufacturer's warranty.

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WARNING: The electrical installations in the operating room in which the unit is connected and operated must comply with the applicable IEC standards.

WARNING: Keep out of reach of patients.

WARNING: The instructions and interface specifications for medical devices and/or system components used in combination must be observed precisely.

WARNING: Combinations of medical devices are only assured to be safe if

- they are identified as such in the respective instruction manuals or
- the intended use and interface specifications of the products used in combination permit this (cf. IEC 60601-1-1 and Chap. 16, IEC 60601-1:2005).

WARNING: Grounding reliability can only be achieved when the equipment is connected to "Hospital Only" or "Hospital Grade" outlet (i.e., approved for use in an operating room environment). Routinely inspect electrical plug and cord. Do not use if inspection reveals damage.

WARNING: To avoid the risk of an electric shock, this unit may only be connected to a power supply network with a protective conductor.

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WARNING: DO NOT use the HAMOU® ENDOMAT® in the presence of flammable anesthetics.

WARNING: Always unplug the unit before all maintenance work (e.g. cleaning).

WARNING: To avoid the risk of electrical shock, do not open the unit. Refer servicing to qualified personnel.

WARNING: Pressure-controlled pump motor. When the operating switch is on, never reach into the pump rollers even if the pump is stationary. The pump may start up at any time when the operating switch is on.

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WARNING: It is important to carefully monitor the fluid level in the irrigation bag to ensure that the irrigant does not completely run out and result in air being pumped into the patient. Another bag of irrigant should be started prior to emptying the bag being used.

WARNING: The HAMOU® ENDOMAT® does not measure intracorporeal pressure. The pressure setting is the pressure at which the irrigant is being pumped into the system. Irrigation line pressure and flow can be read on displays ⑤ and ③, respectively. An intermittent 'beep' tone will be heard if the pressure limit setting is briefly exceeded. For additional safety there is an overpressure valve integrated into the pressure measuring chamber to prevent maximum delivery pressure from being exceeded.

WARNING: Use of electromechanical medical devices can cause electromagnetic interference in other devices, particularly cardiac pacemakers. Precautions must be taken to ensure that the patient's well being is maintained in the event of such interference.

WARNING: DO NOT reuse the disposable tubing set. DISCARD after use.

WARNING: Reusable tubing must be cleaned, disinfected and sterilized before use and subsequent reuse using validated preparation procedures.

WARNING: For all cleaning and disinfection procedures on contaminated instruments, the guidelines of the Employer's Liability Insurance Association and equivalent organizations must be observed.

WARNING: The tubing set must be vented before every examination to remove any air bubbles.

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CAUTION: When replacing fuses, use only fuses of the correct rating.

CAUTION: The HAMOU® ENDOMAT® may only be operated with tubing sets and accessories which have been designated as suitable for the unit by KARL STORZ.

CAUTION: Avoid allowing fluids to enter the unit. Do not store liquids on or above the unit.

CAUTION: Always position the HAMOU® ENDOMAT® at or above the position of the patient to assure proper function of the unit.

CAUTION: If a power failure occurs after the tubing set has been installed on the unit, the self-test procedure can be by-passed by depressing the footswitch when turning on the power switch.

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NOTE: This unit has been marked in accordance with the European Directive on waste electrical and electronic equipment (WEEE). At the end of its useful operating life, dispose of the unit as electronic scrap.

Please ask either KARL STORZ GmbH & Co. KG, a KARL STORZ subsidiary or your specialist dealer for information on your local collection point. Within the scope of application of this Directive, KARL STORZ GmbH & Co. KG is responsible for the proper disposal of this unit.

NOTE: Make certain to observe national regulations/laws regarding disposal (e.g. used irrigation fluid, tubing sets, irrigation bags, etc.).

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Before using the unit on the patient it is imperative that you be acquainted with how the unit operates and is controlled.

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6.2 Intended use

Suction/irrigation pumps and their accessories are used for the introduction of irrigation fluids into organs and operating fields as well as the suctioning off of irrigation and bodily fluids, secretions and tissue during endoscopic interventions.

Use of the unit in fields other than those indicated above is not allowed for safety reasons.

Contraindication

Suction/irrigation pumps without a pressure monitoring system must not be used for the dilation of hollow organs and joints. Serious injuries to the patient may result if this use exclusion is not observed.

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The HAMOU® ENDOMAT® may only be operated with tubing sets and accessories which have been designated as suitable for the unit by KARL STORZ.

No responsibility for the safe operation of the instrument can be accepted if any tube systems other than those specified are used. The following serious defects may arise in particular:

- Inadequate grip on the tube connectors means a loss of sterility and may lead to the unit being damaged by ingress of liquid.
- Premature breakage of the pump tube segment and consequent loss of sterility and damage to the unit caused by ingress of liquid.
- Display of an incorrect actual irrigation rate due to deviations in tube diameter.

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Unauthorized conversions or modifications to the unit are not allowed for safety reasons.

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6.3 User qualification

The HAMOU® ENDOMAT® may only be used by physicians and medical assistants who have a corresponding specialized qualification and who have been instructed in use of the unit.

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6.4 Training in the operation and function of the device

Training should only be provided by persons equipped to do so on the basis of their skills. Your local representative or responsible KARL STORZ member of staff is available to provide training on your premises and to give you information about further training options.

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6.5 Patient profile

6.5.1 Laparoscopy

Gender: irrelevant
Age: irrelevant
Weight: irrelevant
Nationality: irrelevant
State of health: suitable for treatment in the physician's judgment

6.5.2 Gynecology

Gender: female
Age: irrelevant
Weight: irrelevant
Nationality: irrelevant
State of health: suitable for treatment in the physician's judgment

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6.6 Intended region of treatment on the body

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6.7 User profile of physician and assistants

Recognized medical skills in the field of application (specialist physician, qualified medical staff).
Adequate powers of comprehension to rationally assess the situation in hand.
Familiarity with the use of technical apparatus.
Suitable command of one of the languages used by the device and in the instruction manual.
At least one comprehensive session of training in the use of the device.
No physical impairments which limit or prevent perception of acoustic or visual alarm signals.

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6.8 Intended conditions of use

6.8.1 Use

The device is intended for use in hospitals and doctors' offices. The technical data and ambient conditions are described in the instruction manual.

6.8.2 Other intended conditions

Frequency of use: one or more times a day.
Length of use: from a few minutes to several hours a day.
Place of installation: Positioning on a level, vibration-free surface.
Mobility: can be moved if positioned on a cart.
Combination: can be used on the patient at the same time as other equipment required for the operations.
Control: can be controlled via the KARL STORZ SCB.

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
6.9 Safety precautions at the site of installation

The unit may only be used in medical rooms in which the electrical equipment has been installed in accordance with applicable national regulations.

It is not intended for use in hazardous zones. This means, for example, that when using easily combustible and explosive inhalation anesthetics or mixtures thereof, the unit must not be operated inside the hazard zone shown in the diagram. This also applies to highly combustible and potentially explosive chemicals, e.g. skin disinfectants and fast-acting surface disinfectants.

The unit is equipped with a connector for attaching a ground line. It should be connected up in accordance with national regulations.

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6.10 Safety precautions when operating the unit

It is the user's responsibility to make sure the unit is safe and operates properly before using it.

During treatment with the ENDOMAT® the patient must be treated and kept under observation with the usual medical care. This includes keeping a check on the progress of treatment, as well as monitoring the vital levels and the anesthetic effects.

Any treatment may only be performed if there is visual observation of the unit in operation.

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6.11 Safety features

The HAMOU® ENDOMAT® is equipped with the following safety features:

6.11.1 Self-Check

A self-check procedure is performed each time the unit is switched on. The components checked are as follows:

- internal self-check of the electronic circuitry
- tube detection test
- pressure measuring chamber detection test
- pressure sensor test

If the results of any of these tests are negative, it will not be possible to start the unit. In this case the unit will issue an intermittent warning sound for 10 seconds. The error display will continue to flash until the error is remedied.

i NOTE: In order to insert the tube set before putting the unit into operation, the self-check can be bypassed by pressing the footswitch during power up.

Make sure the unit performs a self-check procedure at least once every day of operation.

6.11.2 Display element test

On switching on the instrument all the display elements will light up. This is to enable the user to check the display for any faulty segments.

6.11.3 Monitoring the irrigation pressure

Irrigation will be stopped as soon as the actual pressure value exceeds the specified setpoint.

If the deviation lasts for more than several seconds, there will be an audible warning signal.

As an additional safety measure the pressure measuring chamber was fitted with an integrated overpressure valve which prevents the maximum irrigation pressure from being exceeded.

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
7. Installation and operating instructions

7.1 Unpacking the equipment

Carefully unpack the HAMOU® ENDOMAT® and remove the unit and its accessories from their packaging. Check for missing items and evidence of shipping damage. File any complaints with the manufacturer or supplier immediately.

If possible, retain the original packing materials for later use; these can come in handy if the unit has to be transported.

7.2 Basic equipment

- 1 HAMOU® ENDOMAT® 26331020-1
- 1 SCB connecting cord 200901 70, length 100 cm
- 1 Power cord 400 A
- 3 HYS disposable hysteroscopy tubing kit 031117-01
- 3 LAP disposable laparoscopy tubing kit 031118-01
- 1 VACUsafe Promotion pack Suction (2l) 031020-03
- 1 Instruction manual
- 1 Wallet for instruction manual
- 1 Operator's quick reference guide (self-adhesive)
- 8 Label strips

7.3 Installation

NOTE: The HAMOU® ENDOMAT® including the accessories connected may only be used in medical rooms in which the electrical equipment has been installed in accordance with applicable national regulations.

WARNING: The unit is not intended for use in hazardous zones. Do not operate the unit within the hazard zone shown in the diagram while explosive anesthetic gases are in use.

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Fitting a label strip to the front panel

Undo the left-hand front side cover with a screwdriver.

! NOTE: The appropriate controls/indicators on the unit are designated on the label strip.

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Remove any unmarked label strip. Select a label strip ⑩ (currently 8 languages available) and slide it between the plastic and the metal front.

Fold the strip along the broken line to an angle of 90° and approximately 2 mm further on (toward the end of the strip) fold it back to an angle of 90° so that the side cover can be refitted without any obstructions.

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The unit is equipped with a connector ⑪ for attaching a ground line.

The unit's ground line should be installed by a qualified electrician.

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WARNING: Operate the unit only with the system voltage stated on the rating plate.

Connect power cord. Insert power cord into power cord receptacle 23 as far as it will go.

WARNING: Only insert the power plug into and remove it from electrical outlets located outside areas subject to explosion hazards.

i NOTE: Please ensure that access to the power plug is ensured at the site of installation.

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7.4 KARL STORZ SCB

NOTE: To avoid mistakenly pulling out the SCB connecting cable, the SCB connector possesses a protection device.

Pull back the protection device of the SCB connector and insert the connector into one of the SCB sockets 20. Connect the other end of the cable to the KARL STORZ SCB (KARL STORZ Communication Bus) control unit or other SCB units (see KARL STORZ SCB or SCB control NEO System Instruction Manual).

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If necessary, connect footswitch (see Accessories) to the socket 19.

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7.5 Mode of operation selection

The HAMOU® ENDOMAT® unit is suitable for both hysteroscopic and laparoscopic applications. The appropriate mode of operation will be automatically selected when the correct type of tubing for the application is inserted. The unit determines the intended application from the differing 'encodings' of tubing sets and then automatically selects the appropriate mode of operation.

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7.6 Operating the unit

Switch on the unit at its main power switch ①. The red pilot lamp ⑫ will illuminate. The unit will then conduct a self-check procedure, checking all system components for proper operation.

The indicators on the setpoint bargraphs ④, ⑥ and ⑧ for flow, irrigation pressure and suction pressure are flashing.

i NOTE: Once switched on, the unit awaits confirmation or a change of the last set values used (see note, p. 17). This is signalled by flashing of the corresponding LEDs in the bargraph displays. The LED in the irrigation START/STOP button ⑬ also flashes. As long as the displays are flashing, the START function is blocked.

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7.7 Set value confirmation

To confirm all set values, press one of any ± buttons or the START/STOP button ⑬.

The LED in the START/STOP button ⑬ goes out and the LEDs in the bargraph displays are continuously illuminated.

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7.8 Power-up with irrigation START/STOP button pressed down

At device power-up with the footswitch pressed down or with the irrigation START/STOP button ⑬ pressed down, tests for the following routines are not performed: pressure measuring chamber detection test, irrigation sensor test and irrigation sensor calibration.

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7.9 Inserting the tubes

NOTE: The tubing set must not be inserted until the unit has been switched on with the power switch ①. If the tubing set is inserted before switching on, the self-check procedure blocks the equipment. If the footswitch or the irrigation START/STOP button ⑬ is pressed simultaneously while switching on the power switch, the self-check procedure may be avoided. Make sure the unit performs a self-check procedure at least once every day of operation.

When inserting the tubing set, push button ⑨ to release the pump rollers. The pump rollers spring back.

CAUTION: Use only tubing sets designed for the HAMOU® ENDOMAT®.

NOTE: The pump tube connectors have different widths so that the pump tube can be inserted in one direction only.

7.9.1 Disposable tubing set

Place a length of tubing on the tube holder.

Position the pressure measuring chamber on the pressure transducer ⑪ and turn clockwise until the pressure measuring chamber locks into place (bayonet-lock).

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The instruction manual delivered with the tubing sets must be observed.

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Connect the correct end of the length of tubing to the irrigation-liquid bag (penetration fitting), and to the 'LUER-lock' irrigation-line fitting on the instrument to be employed.

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Once the tubing set and pressure measuring chamber have been correctly installed, the pilot lamp ⑫ will switch from red to green.

The mode of operation selected will be indicated by lamp ⑫.

CAUTION: If the pilot lamp ⑫ remains red, either there is a failure somewhere in the unit, or the tubing set has been incorrectly inserted. If the pilot lamp fails to switch from red to green after the tubing set has been removed and reinserted, then the unit will have to be checked by an authorized technician before it can be used.

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7.10 Removing the tubing set

Press pushbutton ⑨ to release pump rollers.

Remove the tubing from the tube holder.

Loosen and detach the pressure measuring chamber from the pressure transducer by turning it counterclockwise.

Install the reusable tubing set for LAP 26330041 or HYS 26330040, available as an optional extra, in the same way as you would the disposable tubing set.

NOTE: For further information on the maintenance/preparation of pressure chambers/tubing sets, please consult the instruction manual delivered with the reusable tubing set.

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7.11 Connecting the tubing set for suction (optional)

Screw the bacterial filter into the vacuum line fitting ⑮ on the suction bottle and attach the suction line to the filter housing.

CAUTION: To prevent contamination of the unit by penetrating suction fluids, a hydrophobic protection filter must be inserted between the suction bottle and the unit.

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Connect the other end of the suction line to the fitting ⑭.

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Use a length of tubing to interconnect the suction line fitting of the fluid collection bottle ⑯ and the suction line fitting on the instrument to be employed.

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i NOTE: Pressure, flow rate, and suction pressure settings are retained in the unit's memory while the unit is switched off, and will be redisplayed flashing when the unit is switched on. To confirm the set values, one-time activation of any ± buttons or the START/STOP button ⑬ is enough, or reset them as described following.

Then the LED in the START/STOP button ⑬ goes out and the LEDs in the bargraph displays are continuously illuminated.

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7.12 Presetting irrigation parameters

Enter the flow rate setpoint using the associated ± buttons ⑮.

The value set will be indicated on the associated display ④.

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
Enter the limiting irrigation pressure using the associated ± buttons ⑰.

The value set will be indicated on the associated bargraph display ⑥.

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7.13 Presetting suction pressure

Enter the suction pressure setpoint using the associated ± buttons ⑱. The value set will be indicated on the associated bargraph display ⑧.

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7.14 Test for proper operation and air-relief of the tubing system

WARNING: Test this equipment prior to each surgical procedure to ensure that it functions correctly.

NOTE: Have a collecting container ready to catch irrigation liquid as it runs out.

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
Press the irrigation START/STOP button (13) or the footswitch until all the air has been let out of the tubing system (i.e. no bubbles) (see Warning, page 6).

The flow rate indicated on bargraph display (3) must agree with the flow-rate setpoint indicated on display (4).

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Gradually close the metering valve on the instrument. As soon as the pressure indicated on bargraph display (5) rises above the limiting value indicated on display (6), the reel pump should be halted. If the actual value exceeds the set value for a number of seconds there is an audible warning signal.

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Pinch off the suction line with your fingers and check whether suction line underpressure is formed. The suction line pressure indicated on display (7) must agree with the suction line pressure setpoint indicated on display (8).

WARNING: If there are any substantial deviations from these results ($> \pm 10\%$), the unit must be inspected by an authorized technician before further use.

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7.15 Initiating irrigation

Irrigation can now be initiated by actuating the irrigation START/STOP button ⑬, or by depressing the footswitch.

Irrigation line pressure and flow rate can be read on displays ⑤ and ③, respectively.

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If the set pressure should be briefly exceeded the pump will be halted and an audible alarm will be heard after a few seconds.

For passive safety there is an overpressure valve integrated into the pressure measuring chamber to prevent maximum irrigation pressure from being exceeded.

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7.16 Initiating suction

Suction can be initiated by actuating the suction START/STOP button ⑮.

Suction pressure can be read on display ⑦.

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8. Maintenance


8.1 Fuse replacement

Switch off the unit and remove the power plug from the electrical outlet.

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Remove the power fuse holder ② with a screwdriver or other tool.

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! **CAUTION:** Only use fuses of the correct rating.

Insert new fuses of the appropriate rating.

	100...240 V~
Power fuse	2 x T2AL250V

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Refit the power fuse holder ②.
Reconnect the power cord.
Test the unit for proper operation.

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8.2 Reprocessing

8.2.1 Wipe-down disinfection of device

WARNING: Always disconnect from the power supply before cleaning!

CAUTION: Make absolutely sure that no liquid can get into the device housing.

The exterior surfaces of the medical device must be wiped clean with a disposable cloth moistened with disinfectant or a ready to use soaked disinfectant cloth. Due to their protein-fixating effect and material incompatibility, alcohol-based agents must not be used. At the end of the necessary exposure time for the disinfectant as specified by the manufacturer, wipe the surface with a dry and low lint cloth.

8.2.2 General warnings for reusable medical devices

WARNING: Risk of infection: Incorrectly reprocessed medical devices expose patients, users and third parties to a risk of infection as well as the risk that the medical device may malfunction. Observe the manual "Cleaning, Disinfection, Care and Sterilization of KARL STORZ Instruments" and the accompanying documentation.

WARNING: Risk of infection: These medical devices are not sterile when delivered. The use of non-sterile medical devices poses a risk of infection for patients, users and third parties. Inspect medical devices for visible contamination. Visible contamination is an indication that reprocessing has not been carried out, or has been carried out incorrectly. Reprocess the medical devices before initial use and before and after every subsequent use using validated procedures.

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WARNING: When carrying out any work on contaminated medical devices, the guidelines of the Employers' Liability Insurance Association and equivalent organizations striving to ensure personal safety must be observed.

CAUTION: When preparing and using solutions, follow the chemical manufacturer's instructions paying close attention to proper concentration, exposure time, and service life. Prolonged immersion and incorrect concentration may result in damage. Bear in mind the microbiological range of action of the chemicals used.

CAUTION: Danger of damage to medical devices: Medical devices may be damaged by using chemicals that have not been approved by KARL STORZ. Only use chemicals approved by KARL STORZ for preparation. A complete list can be found on the Internet at www.karlstorz.com.

CAUTION: National laws and regulations must be observed.

NOTE: The manual "Cleaning, Disinfection, Care and Sterilization of KARL STORZ Instruments" can be downloaded or requested by visiting www.karlstorz.com.

Reusable tubing

WARNING: Disposable tubing sets, locks and hydrophobic bacteria filters must not be reused.

8.2.3 Accessories for preparation

Accessories required for carrying out preparation:

Brushes: Art. 27652 / 27651AL

Miscellaneous: Art. 27660

8.2.4 Disassembly

The medical device must be disassembled into its individual components as far as possible or opened.

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8.2.5 Preparation for cleaning and disinfection

Heavy soiling, corrosive solutions and pharmaceuticals must be removed from the medical device immediately after use. To this end, preclean the medical device by wiping down and rinsing, for example. As a general rule, KARL STORZ recommends manual precleaning under cold running water.

8.2.6 Manual precleaning Brushing the surfaces

Any heavy soiling of the surfaces must be completely removed by cleaning them mechanically under cold running water with the aid of a brush or sponge.

Brushing the lumina

In order to remove all visible contamination, the working channel, lumina and hollow spaces of the medical device must be precleaned under cold running water with the aid of a suitable brush.

Rinsing out with a pressure water gun

The working channel, lumina and hollow spaces of the medical device must be rinsed using a pressure water gun and a suitable rinsing attachment. Use 3 x 5 second blasts from a pulsed jet of water with a pressure of at least 2.5 bar for rinsing.

8.2.7 Manual cleaning

The medical device must be completely immersed in a cleaning solution. Through disassembling and opening the instrument, it must be ensured that even surfaces with restricted access are covered and that no air bubbles are present. After the necessary exposure time, clean the instrument mechanically with a brush or sponge. Finally, it must be rinsed with cold water to ensure neutralization.

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8.2.8 Manual disinfection

The medical device must be completely immersed in a disinfectant solution. Through disassembling and opening the instrument, it must be ensured that even surfaces with restricted access are covered and that no air bubbles are present. Following the exposure time, the medical device must be rinsed several times with completely demineralized water or microbiologically pure, sterile water in order to remove all chemical residues. Finally, all of the surfaces, channels, and lumina are dried completely with sterile compressed air. The cleaning gun with accessories (Art. No. 27660) is ideal for this purpose.

Suction bottle and cap

Carefully pull the tubes from the suction lances and take the bottle out of the holder.

i NOTE: Prompt transportation to the preparation station and a subsequent decontamination process must be ensured.

Preparation station / CSSD: Carefully remove the cap from the suction bottle, empty the contents into the sink and rinse out, repeat rinsing if necessary.

Rinse out the suction bottle with a disinfectant solution and then fill it with the remaining disinfectant solution. Immerse the cap in a plastic basin with disinfectant solution making sure that all air bubbles have been removed and, if necessary, rinse out/fill the attachments/connectors with disinfectant solution using an injection aid.

For further procedures see section 8.2.8

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8.2.9 Machine cleaning and disinfection

The following methods for machine decontamination have been validated and approved subject to compliance with the process parameters described in the manual "Cleaning, Disinfection, Care and Sterilization of KARL STORZ Instruments":

- Machine cleaning
- Thermal disinfection

Thermal disinfection is preferable. The relevant national requirements and the A_0 value must be taken into account when using this method.

A suitable slide-in tray or instrument holder, to ensure that the medical device is thoroughly rinsed out or through, must be selected following consultation with the manufacturer of the device.

i NOTE: If necessary, the instrument must be dried off afterwards by hand.

8.2.10 Connecting

In order to ensure effective machine cleaning and disinfection, the instrument must be connected up to the washer and disinfectant. It must be ensured that the instrument is rinsed out by means of appropriate rinsing lances

Suction Bottle

i NOTE: The bottles must be rinsed out once and precleaned before they can be loaded into the washer and disinfectant.

Carefully place the suction bottle in a suitable tray holder or connect up (to a rinsing jet). Contact between the glass and other medical devices which could result in damage to the material must be prevented during the decontamination program. It is recommended that the suction lances are connected to the rinsing tubes of the washer and disinfectant.

Start the program (93 °C device-specific, alkaline cleaner pH value ≤ 11.5).

Perform a visual inspection of the suction bottle after machine preparation.

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8.2.11 Assembly, inspection and care

The cleaned and disinfected medical device must be visually inspected for cleanliness, completeness, damage and dryness.

- If residues or contamination are still present, the medical device must be manually cleaned and subjected to a full cleaning and disinfection procedure once more.
- Damaged medical devices must be withdrawn from use.
- Afterwards, a functional check must be carried out.

! **NOTE:** During care procedures, use items from the catalog "Care, Sterilization and Storage Techniques".

Reusable tubing

! **WARNING:** Check the entire tubing set for leaks after preparation. This makes it possible to detect a leak early, through which liquid could escape and enter the device.

To check for leaks use leakage tester 13242 XL, for example. To do this, connect one tube end to the leak tester and block off the other end, with your finger, for example. If there is a leak in the tubing set, a pressure drop will be observed on the pressure gauge of the leakage tester.

8.2.12 Packaging systems

Only standardized and approved packaging materials or systems may be used (EN 868 Parts 2-10, DIN EN ISO 11607 Parts 1 + 2, DIN 58953).

8.2.13 Sterilization

The procedures as well as the process-relevant parameters for the individually validated methods are described in detail in the manual "Cleaning, Disinfection, Care and Sterilization of KARL STORZ Instruments". The procedure must be selected taking the respective applicable national requirements into account and in consultation with the device manufacturers.

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The following sterilization methods have been validated and approved by KARL STORZ for this medical device:

Steam sterilization in the fractionated prevacuum procedure

This medical device must be sterilized in its fully assembled state using a fractionated prevacuum procedure (DIN EN ISO 17665-1) at 134-137°C with a minimum exposure time of 4 minutes and maximum of 18 minutes. These methods are only suitable for thermostable instruments.

Suction bottles: Only cleaned, disinfected, dry and functional suction bottles in suitable packaging may be sterilized. The suction bottle and cap can be steam sterilized at 134°C (± 3 °C), hold time at least 4 minutes, fractionated prevacuum method.

i NOTE: To improve sterility efficiency, sterilize the suction bottle in its disassembled state. We recommend sterilizing the suction bottle separately, i.e. packed in sterile goods non-woven material/paper and separated from the other medical devices to be sterilized.

Repreparation limits

The product's service life and correct functioning are largely determined by mechanical stress and chemical influences within the scope of reprocessing and application.

8.2.14 Storage

Store the suction bottles and caps in suitable sterile goods packaging in a closed cupboard to protect them from dust, moisture and temperature fluctuations.

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8.3 Maintenance and safety check

8.3.1 Maintenance

Preventive maintenance is not essential. Regular maintenance can, however, contribute to identifying potential problems before they become serious, thus enhancing the instrument's reliability and extending its useful service life. Maintenance services can be obtained from your local representative or from the manufacturer.

8.3.2 Safety checks/repeat inspections as per IEC 62353

WARNING: Regardless of the accident prevention regulations or testing intervals for medical devices prescribed in different countries, with this device safety checks/repeat inspections as defined by IEC 62353 must be performed and recorded by a qualified electrician at least once a year.

Visual inspection

1. Check the device and accessories for any mechanical damage which may impair functionality.
2. Check that inscriptions relevant to safety are legible.

Electric measurements

- Inspect the device safety fuses
- Measure protective ground resistance in accordance with IEC 62353: The limit values can be found in the current standard.
- Measure earth leakage current in accordance with IEC 62353: The limit values can be found in the current standard.
- Touch current measured according to IEC 62353: The limit values can be found in the current standard.
- Patient leakage current measured according to IEC 62353 (with mains on applied part): The limit values can be found in the current standard.

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Test for proper functioning

Carry out the functional check as per the instruction manual (see section 7.14).

Documentation

The safety check/repeat inspection and results must be documented.

WARNING: If any defects or shortcomings which could endanger patients, personnel or others are ascertained during these safety checks/repeat inspections, the device must not be operated until such time as these defects or shortcomings have been eliminated by qualified technical servicing.

NOTE: For detailed information on the range and implementation of the safety check/repeat inspection, please refer to the latest issue of the relevant Service Manual.

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8.4 Servicing and repair

Defective items of equipment must be serviced and repaired exclusively by persons authorized by us; all repair work must employ original KARL STORZ parts only.

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8.5 Disposal

This unit has been marked in accordance with the European Directive on waste electrical and electronic equipment (WEEE).

At the end of its useful operating life, dispose of the unit as electronic scrap.

Please ask either KARL STORZ GmbH & Co. KG, a KARL STORZ subsidiary or your specialist dealer for information on your local collection point.

Within the scope of application of this Directive, KARL STORZ GmbH & Co. KG is responsible for the proper disposal of this unit.

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8.6 Repair program

For fiberscopes and equipment, individual repair is necessary. Usually to bridge the repair period, you will receive a unit on loan which you then return to KARL STORZ as soon as you receive the repaired unit.

In Germany you can refer repairs directly to

KARL STORZ GmbH & Co. KG
Abt. Reparaturservice
Take-off Gewerbepark 83
78579 Neuhausen

In other countries please contact your local KARL STORZ branch or authorized dealer.

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8.7 Important information

In order to prevent infection, it is strictly forbidden to ship contaminated medical devices. All medical devices must be decontaminated on site, to avoid contact and aerogenous infections (among staff).

We reserve the right to return contaminated instruments/equipment to the sender.

Repairs, changes or expansions which are not carried out by KARL STORZ or by experts authorized by KARL STORZ will invalidate all guarantee rights. KARL STORZ gives no guarantee on the correct functioning of equipment or instruments which have been repaired by unauthorized third parties.

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8.8 Limitation of liability

KARL STORZ shall be liable for the safe operation, operational reliability, and performance of this equipment only subject to the following conditions:

- all assembly operations, system expansions, readjustments, modifications, or repairs have been performed by a person or persons duly authorized by KARL STORZ;
- all electrical installations at the location in which the unit is connected and operated meet the applicable laws and standards; and
- the unit has been used in accordance with its operating instructions at all times.

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8.9 Manufacturer's warranty

The guarantees provided can be found in the Standard Conditions of Business of KARL STORZ.

The medical device must always be sent to your local subsidiary (see "Subsidiaries" section), even during the warranty period.

Opening the equipment or performance of any repairs or modifications to the equipment by unauthorized persons shall relieve us of any liability for its performance. Any such repair or modification performed during the warranty period shall void all warranty.

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9. Technical description

9.1 Troubleshooting

WARNING: Always unplug the unit before all maintenance work.

Symptom:

- Complete failure of the unit.

Possible causes:

- No power from the power line.
- Defective power fuse.
- Power cord is not properly connected to power cord connector.

Remedy:

- Check that there is electricity to the wall outlet.
- Change fuses as described in the instruction manual. Use only fuses of the correct rating.
- Push power cord firmly into power cord connector.

Symptom:

- Poor suction.

Possible causes:

- Leakage occurring in suction lines.

Remedy:

- Check all suction lines, and replace any exhibiting evidence of leakage.
- Check the seating of the sealing cap.

Symptom:

- No suction.

Possible causes:

- Floats of overflow protection unit blocking suction inlet.

Remedy:

- Check the fluid level in the vessel, and empty, if required.
- Clean the floats and check them for freedom of movement.

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9.2 Error codes

If an error occurs when the unit is powered up, the unit will issue an intermittent warning sound for 10 seconds. The error display will continue to flash until the error is remedied. The number of flashing pulses indicates the error code.

Only one error code will be displayed: the one which was detected first. In addition, the error code is shown on display ③ 'Irrigation flow act. value'. The error code corresponds to the number of LEDs lit up. This code is for the user's information or provides Service with preliminary information.

Error	Description of the error codes
1	Pressure measuring chamber attached while switching on
2	Checksum error of the EPROM (checksum ≠ 0)
3	Error during RAM test
4	Error during initial flow measurement, flow value > 0 (≠ 0)
5	Offset error of pressure sensor, p measured < -20 mmHg
6	Offset error of pressure sensor, p measured > 40 mmHg
7	Motor control module too hot > 135 °C
8	Excessive motor current, irrigation pump > 200 mA ± 50%
9	Excessive motor current, suction pump > 200 mA ± 50%
10	Offset error, suction pressure too low > 40 mbar
11	Offset error, suction pressure too high < -40 mbar
12	Tubing detection operated
13	Motor control too hot (during operation)
14	Excessive irrigation-pump-motor current (during operation)
15	Excessive suction-pump-motor current (in operation)
16	Maximum number of SIO errors reached

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- Type of protection against electric shocks:
Protection Class I
- Degree of protection against electric shocks:
Applied part of type BF

According to IEC 60601-1-2:

Observe the information on electromagnetic compatibility in chapter 11 (p. 42-55).

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9.3.2 Directive compliance (for 263310 20-1)

According to Medical Device Directive (MDD):

This medical product belongs to Class IIb

This medical product bears the CE mark in accordance with MDD 93/42/EEC.

i NOTE: The code number after the CE mark indicates the responsible notified body.

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9.4 Technical documentation

On request, the manufacturer will provide those circuit diagrams, itemized parts listings, descriptions, sets of adjustment instructions, and other items of available documentation to suitably qualified user personnel duly authorized by the manufacturer for their use in repairing those components of the unit that have been designated by their respective manufacturers as repairable.

Supply of such technical documentation relating to the unit shall not be construed as constituting manufacturer's authorization of user's personnel, regardless of their level of technical training, to open or repair the unit.

Explicitly exempted herefrom are those maintenance and repair operations described in this manual.

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
We reserve the right to make engineering modifications in the interest of promoting technological progress and generating performance improvements without obligation on our part to submit prior notice thereof.

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10. Spare parts, recommended accessories

10.1 List of Spare parts

Item	Cat. no.
Disposable tubing kit HYS, pack of 10	031117-10*
Disposable tubing kit LAP, pack of 10	031118-10*
SCB connecting cord, length 100 cm	20 0901 70
Power fuses T2AL250V, pack of 10	2027590
Power cord (grounded)	400 A
"Hospital Grade" Power cord (USA)	400 B
Instruction manual	96116008 DER
Operator's quick reference guide	96116008 X

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Postfach 230

78503 TUTTLINGEN

For detailed information please contact:

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11. Electromagnetic Compatibility (EMC) Information

WARNING: Medical electrical equipment needs special precautions regarding Electromagnetic Compatibility (EMC). Observe the EMC instructions in this chapter during installation and operation.

The HAMOU® ENDOMAT® model 26331020-1 corresponds to IEC 60601-1-2:2007 [CISPR 11 Class B] and therefore meets the EMC requirements of the Medical Device Directive (MDD) 93/42/EEC.

These limits are designed to provide reasonable protection against the typical electromagnetic interference to be expected in a medical environment. The HAMOU® ENDOMAT® model 26331020-1 is a Group 1 unit (as per CISPR 11). Group 1 contains all the 'equipment and systems which generate or use RF energy only for their internal functioning'.

NOTE: The tables and guidelines that are included in this chapter provide information to the customer or user that is essential in determining the suitability of the device or system for the electromagnetic environment of use, and in managing the electromagnetic environment of use to permit the device or system to perform its intended use without causing interferences in other devices or systems or non-medical electrical devices. If this equipment does cause harmful interference with other devices, the user is encouraged to try to correct the interference by one or more of the following measures:

- reorient or relocate the receiving device
- increase the separation between the equipment
- connect the equipment into an outlet on a circuit different from that to which the other device(s) is connected.

If you have any further questions, please contact your local representative or our service department.

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Electromagnetic Compatibility (EMC) Information

WARNING: The HAMOU® ENDOMAT® model 26331020-1 should not be used adjacent to or stacked with other equipment. If adjacent or stacked use is necessary, the equipment or system should be observed to verify normal operation in the configuration in which it will be used.

WARNING: The use of portable and mobile RF equipment may have an impact on this and other pieces of medical equipment.

WARNING: The use of accessories and cables other than those specified in the KARL STORZ instruction manual may result in increased emissions or decreased immunity of the HAMOU® ENDOMAT® model 26331020-1. The accessories and cables listed below have been shown to comply with the requirements of IEC 60601-1-2. When using accessories and cables other than those specified here, it is the responsibility of the user to ensure that they comply with IEC 60601-1-2.

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