

KODAK 2200 Intraoral X-ray System

Installation & Service Manual



This document was originally written in English
Revision date: 09/2007

The brand names and logos reproduced in this guide are Copyright.

KODAK is a trademark of KODAK used under License.

RVG, Trophy, are registered trademarks of Carestream Health, Inc. All other names or products referred to in this document are used only for the purpose of identification and may be the trademarks or the registered trademarks of their respective owners.

The RVG technology is the subject of an international patent Registered by Carestream Health, Inc.

Table of Contents

Section 1 : Safety and Regulatory Information

- 1 – Document Conventions**
- 2 – General Safety Guidelines**
- 3 – Warnings and Safety Instructions**
- 4 – Labeling Summary**
- 5 – IEC Symbols Used**
- 6 – Regulatory Information**

Section 2 : System Overview

- 1 – Packaging**
- 2 – Components**
- 3 – Operating Controls**

Section 3 : Installation Layouts

- 1 – Wall Framework Layouts**
- 2 – Optional Configurations**

Section 4 : Checking the Premises Before Installation

- 1 – Equipment and Tools Needed for Installation (Not Supplied)**
- 2 – Required Mechanical Specifications**
- 3 – Required Electrical Specifications**
- 4 – Preparation of the Room**

Section 5 : Installation

- 1 – Preparation for Installation**
- 2 – Mechanical Installation**
- 3 – Timer Configuration**
- 4 – Electrical Installation**

Section 6: Setup

- 1 – Checking the Voltage**
- 2 – Checking the mA**
- 3 – End of Installation**

Section 7: Preventive Maintenance

- 1 – Cleaning and Disinfecting**
- 2 – Maintenance Checks**

Section 8: Corrective Maintenance

- 1 – Main Components**
- 2 – Diagnosis**
- 3 – Part Replacement Procedures**

Appendix

- 1 – Dimensional Diagrams**
- 2 – Technical Specifications**
- 3 – Spare Parts List**

SECTION 1

SAFETY AND REGULATORY INFORMATION

The information contained herein is based on the experience and knowledge relating to the subject matter gained by Carestream Health, Inc. prior to publication. No patent license is granted by this information.

Carestream Health reserves the right to change this information without notice, and makes no warranty, express or implied, with respect to this information. Carestream Health shall not be liable for any loss or damage, including consequential or special damages, resulting from any use of this information, even if loss or damage is caused by Carestream Health negligence or other fault.

1 - Conventions Used in This Manual

CAUTION

Caution points out procedures that you must follow precisely to avoid damage to the system or any of its components, yourself or others, loss of data, or corruption of files in software applications.

NOTE

Notes provide additional information, such as expanded explanations, hints, or reminders.

IMPORTANT

Important highlights critical policy information that affects how you use this manual and this product.

2 - General Safety Guidelines

- This product is designed and manufactured to ensure maximum safety of operation. Operate and maintain it in strict compliance with the safety precautions and operating instructions contained in this manual.
- This product meets all the safety requirements applicable to medical equipment. However, anyone attempting to operate the system must be fully aware of potential safety hazards.
- There are no user serviceable parts in this system. The product must be installed, maintained, and serviced by qualified service personnel according to the procedures and preventive maintenance schedules in the product service manual. If your product does not operate as expected, contact your Service Representative.
- Do not modify this product in whole or in part without prior written approval from Carestream Health.
- The assembly, extensions, adjustments, modifications, and repairs must be performed by an authorized Service Representative. Your radiology system must be installed in premises that comply with applicable standards.
- Personnel operating and maintaining this system should receive training and be familiar with all aspects of operation and maintenance.
- To ensure safety, read all user manuals carefully before using the system and observe all Caution, Important, and Note callouts located throughout the manual.
- Keep this manual with the equipment.

-
- Reading this manual does not qualify you to operate, test, or calibrate this system.
 - Unauthorized personnel are not allowed access to the system.
 - If the product does not operate properly or fails to respond to the controls as described in this manual:
 - Follow the safety precautions as specified in this manual.
 - Stop using the equipment and do not make or authorize any changes to it.
 - Immediately contact your Service Representative, report the problem, and await further instructions.
 - X-ray systems manufactured by Carestream Health comply with safety standards throughout the world for optimum protection against radiation risks.
 - Be aware of the product specifications and of system accuracy and stability limitations. Consider these limitations before making any decision based on quantitative values. If you have any doubts, consult your Sales Representative.

CAUTION

X-rays can be dangerous if used incorrectly. Take precautions even when following the instructions in this manual.

Use conventional, commercially available equipment to protect yourself and your patients against the risk of scattered radiation.

- If you fail to comply with these instructions, Carestream Health will not be responsible for the safety reliability or characteristics of the equipment.

3 - Warnings and Safety Instructions

CAUTION

Do not operate the equipment in the presence of explosive liquids, vapors, or gases. Do not plug in or turn on the system if hazardous substances are detected in the environment. If these substances are detected after the system has been turned on, do not attempt to turn off the unit or unplug it. Evacuate and ventilate the area before turning off the system.

DANGER: THIS IS AN ELECTRICAL UNIT. DO NOT EXPOSE IT TO WATER SPRAY. SUCH ACTION MAY CAUSE AN ELECTRICAL SHOCK OR A MALFUNCTION OF THE UNIT.

WARNING

The user is responsible for the operation and maintenance of this unit.

This unit must be operated only by legally qualified persons.

The cover of the unit must not be opened by the operator.

Inspection and maintenance operations should be carried out only by an approved technician.

WARNING

This unit must be installed in an x-ray room that complies with current installation standards. From this location, visual or audio communication must be maintained with the patient, together with access to the control interface during exposure.

WARNING

Do not operate the unit if there is the threat of an earthquake.

Following an earthquake, ensure that the unit is operating properly before using it again.

Failure to observe this precaution may expose patients to hazards.

WARNING

X-ray equipment can be hazardous to patients and the operator if the exposure safety factors and operating instructions are not observed.

WARNING

Do not place objects within the field of operation of the unit.

WARNING

We recommend that the patient and the operator wear protective lead-lined aprons, unless other Radiation Protection Protocols apply locally.

Ensure that any parts of the unit that may come into contact with the patient and the operator have been disinfected after each patient has been exposed to x-rays.

If the unit develops a fault, turn it off (O) and display a sign that states "Out of Service."

WARNING

The operator must ask the patient to refrain from moving during the entire period of exposure.




In the European Union, this symbol indicates that when the last user wishes to discard this product, it must be sent to an appropriate facility for recovery and recycling.

Contact your local sales representative for additional information on the collection and recovery programs available for this product.

4 - Labeling Summary

Safety Labels



CHASSIS GROUND STUD



ATTENTION: CONSULT ACCOMPANYING DOCUMENTS



CAUTION: IONIZING RADIATION

5 - IEC Symbols Used

The system may have labels with one or more of the following symbols. These symbols indicate the IEC standards to which the system conforms.



Caution — consult accompanying documents



Protective earth



Power ON



Power OFF

6 – Regulatory Information

The product conforms to the following safety standards:

- IEC/EN 60 601-1 Medical Electrical Equipment General Requirements for Safety,
- IEC/EN 60 601-2 Medical Electrical Equipment Electro-Magnetic Compatibility Requirements and Tests.

CE Conformity

This product conforms to the requirements of EU Council Directive 93/42/EEC. The Kodak Intraoral X-ray System is a Class II b medical device, which bears the following mark of conformity:



U.S. Regulations


CAUTION

U.S. federal law restricts this device to sale by or on the order of a dentist.

SECTION 2

SYSTEM OVERVIEW

1 - Packaging




Box dimensions (L x W x H):

1 000 x 600 x 300 mm (39-3/8 x 23-5/8 x 11-13/16 in.)

The Kodak 2200 unit is composed of the following components arranged in special expanded polystyrene packing to guarantee complete protection during transport, as shown in the diagram above.

N°	DESCRIPTION
1	Complete generator
2	Scissor arm with generator cable
3	Plastic covers for the scissor arm
4	Extension arm length: 470 mm (18-1/2 in.), 648 mm (25-1/2 in.) or 825 mm (32-1/2 in.) as specified in the order
5	Wall framework
6	Control timer
7	Power board
8	Accessories: brake for pivot, brake for arm...
-	Technical documentation

2 - Components



Your intraoral x-ray system unit is composed of:

- A high frequency x-ray generator comprising:

- Ø a transformer and associated electronics, and an oil-bathed x-ray emitting tube

- Ø a beam-limiting device that limits radiation to a diameter of 6 cm (2-3/8 in.) on the skin and ensures a distance of 20 cm (7 7/8 in.) between the skin and the x-ray tube focal spot

- Ø an angulations scale and a handle to facilitate positioning


- , A wall framework containing the high frequency generator's control power board, designed to support its mechanical stand

- f* A timer/control unit for the x-ray generator, with the following features:

- anatomical selection and digital display of parameters (kV, mA, exposure time)
- a self-test of the microprocessor each time the unit is activated
- an alarm in the event of incorrect operation
- a RVG key that automatically adjusts the exposure parameters (time and mA) if you are using a Kodak RVG (RadioVisioGraphy) sensor

„ A scissor arm:

The scissor arm makes it possible to position the generator easily and precisely. It is wall-mounted with an extension of 470 mm (18-1/2 in.), 648 mm (25-1/2 in.) or 825 mm (32-1/2 in.).



EXTENSION	R	SPAN A
CG645	470 mm (18-1/2 in.)	1 700 mm (66-15/16 in.)
CG646	648 mm (25-1/2 in.)	1 880 mm (74 in.)
CG648	825 mm (32-1/2 in.)	2 050 mm (80-11/16 in.)


... On/off switch with built-in light.

† rectangular collimator of different sizes adapted to films and RVG sensors (optional).

Options

- Separate control timer
- Separate exposure switch if the control panel is mounted on the wall
- Ceiling-mounted unit
- Mobile stand
- Floor-column.

3 - Operating Controls




- 1 LCD display
- 2 kV selection
- 3 RVG function
- 4 X-ray emission indicator light
- 5 Radiography control
- 6 Receptor sensitivity selection and adjustment
- 7 Adult / Child selection
- 8 Maxillary programming
- 9 Mandible programming
- 10 Occlusal function
- 11 Bitewing functions

SECTION 3

INSTALLATION LAYOUTS

1 – Wall Framework Layouts

There are 3 different ways to position the wall framework: Right, Left or Vertical position. These 3 different layouts can apply to the standard mounting (as shown below), or with the optional configurations with separate timers. In these cases, framework and control timer must be separated by 15 meters (16.4 yards) maximum.



2 – Optional Configurations




CAUTION

To be achieved, all these optional configurations require the addition of the remote timer box (**Cat 5154315**) or the separated exposure switch (**Cat 5154307**) to the standard configuration. These accessories have to be ordered separately. They do not come with the standard kit.

2-1 – Separate control timer on the wall

In this configuration, the timer is mounted on the wall on a dedicated wall-mount assembly, at a distance from the framework up to 15 meters (16.4 yd.).



2.2 – Separate exposure hand-switch




1 – connection box

2 – exposure switch

NOTE

In this configuration, the exposure is activated only from the exposure switch. The exposure button on the control timer is disabled.

2.3 – Dual-command switches with separate control timer and exposure switch (mainly for U.S. market)



1 – timer = 1st button

2 – switch = 2nd button

NOTE

In this configuration, the exposure is activated by acting simultaneously on the exposure switch and on the exposure button from the control timer.

SECTION 4

CHECKING THE PREMISES BEFORE INSTALLATION

1 - Equipment and Tools Needed for Installation (Not Supplied)

- 1 Multimeter with range 300 V-AC ~ 1% and 30 V-DC = 1%, internal resistance greater than 100 k Ω , equipped with 2 measurement wire grips
- 1 Tape measure
- 1 Socket wrench of 7 mm (0.28 in.)
- 1 Set of metric Allen keys
- 1 TorX key reference # 10
- 1 Spirit level
- 1 Plastic mallet
- 1 Impact drill with a set of drill bits from 3 to 13 mm (.12 to .51 in.) suitable for the material of the wall
- 1 Small screwdriver (3 mm) (.12 in.) for electrical connections
- 1 Medium screwdriver
- Some plastic clamps for attaching cables
- Three-wire power supply cable (2 conductors + ground). Recommended cross-section 2.5 mm² (14 AWG), length not exceeding 12 meters (13 yd.) in 100-130 V and 24 meters (26 yd.) in 230-240 V. For greater lengths, use a cable with a bigger cross-section corresponding to regulations in force in the country concerned.
- Shielded cable 2 x 0.22 mm² (2 x 24 AWG) for the separate exposure switch (optional)

To mount the wall framework, use mounting tools suitable for the type of wall concerned (see "Required Mechanical Specifications" below).

2 - Required Mechanical Specifications



WARNING

The Kodak 2200 unit must be installed so that it is impossible to rotate the scissor arm through 360°.

Otherwise, the power supply cable can fail and cause electrical damage.

The wall framework which carries the x-ray unit must be strongly attached to the wall. Choose an attachment system suitable for the type of wall and capable of withstanding a tear-off force of **147 kg (324 lbs)** per anchor point.

Check that the wall is flat and level to avoid any drift in the arm position. Use spacers if the wall surface is not level.

3 - Required Electrical Specifications

A dedicated three-wire electrical line protected by a 16 A breaker (D curve) and a RCCB of 30 mA from the electrical panel must be provided.

A power supply line generally consists of a three-wire cable (2 conductors + ground). The minimum cross-section is 2.5 mm² (14 AWG) and the wire colors and characteristics correspond to each country's electricity regulations. The maximum lengths are 12 m (13 yd.) in 100-130 V and 24 m (26 yd.) in 230-240 V. For greater lengths the cross-section will be increased proportionally, e.g., 3 mm² (12 AWG) for a length of 32 m (35 yd.) in 240 V.

The line's apparent resistance must be less than or equal to 0.2 Ω in 110-130 V and 0.5 Ω in 230-240 V.

The x-ray unit must have a fixed connection to the electrical power supply network.

Rated voltage (off-load)	Minimum	Maximum	Maximum line current
100 V, 110 V, 130 V	90 V	144 V	12 A
230 V, 240 V	207 V	264 V	5 A

IMPORTANT

The maximum line current is obtained by following the instructions for connecting the generator's power supply cable (see paragraph 4.2.2 in Section 5).
The maximum on-load voltage variation on the line must not exceed 3%.



WARNING

If other appliances are installed on the same line, interference and voltage variations may cause your x-ray unit to function abnormally.
To prevent electric shock, this equipment must be ground connected.

Preparing the Room

Compliance with all national and local codes, as well as manufacturer's specification (given above in Sections 3 and 4), is mandatory for high and low voltage wiring.

See Appendix 1 for dimensional diagrams.

SECTION 5

INSTALLATION

IMPORTANT

The unit is supplied with the wall framework oriented to the right for a standard installation.

1 - Preparation for Installation

1.1 - Preliminaries

Only qualified technicians are authorized to install the Kodak 2200 unit, and it must be installed in compliance with the mechanical and electrical installation instructions defined below.

Regardless of the type of installation, the timer must be installed so that the patient and the selected parameters can be seen from the point of operation, and at a minimum distance greater than 2.5 meters.

CAUTION

Make sure that the various cables required for installation have been put in place. Check that the voltage shown on the label of your generator corresponds to the voltage provided by your electrical power supply.


Check that the position of the jumpers on the wall framework's power board corresponds to the voltage, see paragraph 3.1 in this Section.




1.2 - Framework orientation downward or to the left

Unpack and remove the plastic cover from the framework.

Remove the pivot arm support from its initial position using a 6 mm Allen key to remove the two Socket Head Cap screws (a) and washers serving as an axis, and then the two fixing screws (b).



Remove the support clamp (c) using a 4 mm Allen key and take the clamp out through the back of the framework.



Place the clamp, passing it through the back of the framework in the openings corresponding to the orientation of the wall framework: horizontal framework oriented to the left, or framework oriented downward.

Fix the clamp with its screw. Tighten firmly.

Position the pivot arm support in the clamp so that it is held securely by the guiding ribs.

Fix the arm support with its axis screws (a) and washers and the other 2 screws (b).

Tighten gently.

2 - Mechanical Installation

2.1 - Wall framework

Remove the plastic insulator from the chassis, and save it for later (see paragraph 4.2.1).

It is advisable to install the top of the framework approximately 1.20 m (47-1/4 in.) above the floor.

A pre-installation template for easier positioning is available on the top flaps of the unit packaging. If you choose to use it, cut the cardboard properly.

Position the wall framework or the template at the required point with its chosen orientation (see Section 1, paragraph 1.2). Use a spirit level to check that it is horizontal (or vertical in vertical configuration).


Mark the 4 mounting points on the wall.

Drill the holes in the wall and install the fixing system chosen according to the type of wall.

Pass the power supply cable through the holes in the framework ((a) if cables are in the wall, (b) if cables are on the wall), as well as the RVG link cable and the control cable for the remote timer if necessary.

Mount the framework to the wall (c), and if necessary insert shims to ensure levelling of the frame.

Screw and tighten firmly, having first checked that the wall framework is perfectly vertical or horizontal.



2.2 - Extension arm


Check that each pivot of the extension arm is fitted with a thin nylon ring (a).

DO NOT LUBRICATE THE PIVOTS. The scissor arm and the wall framework are fitted with self-lubricating rings.



WARNING

Install the extension arm correctly, with the stop at the scissor arm end (c).



Mount the extension arm's pivot in the top hole of the wall framework.

Install the pivot's brake kit (b).

The parts in this kit must be assembled in the correct order to ensure optimum operation of this system.

Tighten the screw only moderately at first.


2.3 - Scissor arm



WARNING

Do not loosen the strap holding the scissor arm until the generator has been mounted to it. Otherwise, if the arm suddenly swings back it could be damaged and might injure the operator. The scissor arm has not been designed for use without a generator.

The scissor arm ring is self-lubricating, and so **there is no need to lubricate** the extension arm pivot.




- Take the generator's power supply cable coming out of the scissor arm's axis and feed it into the extension arm pivot, bringing it out underneath.

, Install the scissor arm on the extension arm while sliding the power supply cable.

f Pass the cable along inside the extension arm using the opening located under the arm and bring it out on the wall side, through the opening located on the top of the extension arm.

„ Slide the cable into the extension arm pivot and bring it out in the wall framework.

... Make a loop with the cable and attach it to the pivot bloc with the cable clamp, leaving enough length to connect. Slide the extra length inside the extension arm.



† Install the two covers for the power supply cable on the extension arm, one on the wall framework side, and the other on the scissor arm side.

2.4 - Generator


- You will find the generator's blocking key (A) in the accessory box.

, Make sure the ring preventing full rotation (B) is in place.
Remove all excess grease from the tube with a paper towel prior to putting the connector in place.

f Bring the generator (C) from underneath, with the power supply socket toward the top (D).

„ Push upwards until the scissor arm and the generator are in contact.

... Keep the generator in this position.
Position the blocking key in its slot.
You can now remove your hand from the generator.



H214_0002GC

† Insert the plug in the socket located inside the scissor arm and push until they lock together.

‡ Install the two plastic covers (E) and (F) on each side of the arm. Make sure they are properly positioned.


^ Remove the strap holding the arm.

2.5 - Installing the brake

The brake prevents drift when the arm is extended.

- Open the scissor arm slightly to install its brake kit (A).

, The parts in this kit must be assembled in the correct order to ensure optimum operation of this system. Tighten the screw only moderately at first.




H214_0003AC



2.6 - Mechanical adjustments

2.6.1 - Adjusting the arm's horizontal positioning

Extend the arm fully, perpendicular to the wall framework, and place a spirit level as shown in the figure below.



H214_0004HC

Adjust the stop screw • so that the extension arm is horizontal. Maintain the extremity of the arm while operating to avoid constraints. Check that whatever position the arm is in, it does not drift.

If necessary slightly loosen the 2 screws working as an axis , and the 2 screws f in order to facilitate adjustment.

Tighten the 2 screws f moderately, but firmly tighten the axis screws , .

2.6.2 - Adjusting the arm drift

Remove the brake cover on the scissor arm side.

Adjust the brake to prevent horizontal drift of the arm in all positions, while conserving flexibility of movement.


Screw moderately to prevent drift.

Unscrew to make the movement more flexible.

This adjustment should be carried out on both brakes, in the wall framework and on the scissor arm.

WARNING


Never tighten the brake screws firmly as this will destroy the brake pads.



H214_0005AC

Put the cover back on the scissor arm side.

5-9





Power board layout	
Position	Description
1 (or B)	Configuration 110V / 230 V (2 different locations on the power board))
2 (or A)	ST1 configuration (used only for 110V)
3	12V led D48 and PT19 - 5V led D51 and PT20
4	Dipswitches SW2 (by default, all switches are set on OFF)
5	KV return PT13 - IHEAT return PT15
6	Neutral fuse (except for mobile option)
7	“Ready state” light connector 230 V – 60 W mx
8	Fuse: 5 A for 230 V - 10 A for 110 V
9	Main power supply connector
10	Separate x-ray switch connector
11	Synchro RVG connector
12	I2C bus PT22 PT23 PT24
13	Vdac ref. PT21
14	KV ref. PT14
15	12V non isolated: D65 and PT10
16	H Bridge control PT11 and PT12
17	VHEAT PT16 and D41
18	MA return PT17 Relay R control PT18
19	Arm cable connector
20	J1 control timer connector
21	Ground - 0 mA PT25
22	Connection between main power board and adapter board for the control timer cable


3.1 - Power board configuration

Check the main voltage through the correct position of the jumpers on the board in the 2 positions (A) and the position (B).

Check that the sets of jumpers on the board correspond to the measured voltage.

- on 110 V: 2 times 2 jumpers (B)
1 time 1 jumper (A)

- on 230 V: 2 times 1 jumper (B)
no jumper (A)



Verify that fuse F2 (6) is a neutral link, except for the optional mobile stand unit.

3.2 - Power board installation


3.2.1 - Mechanical installation

Check that the electrical cable is not live.

Position the various cables in the framework using the cable-clamps (b) and cable ties (a) (not supplied), according to the diagram below:

- Generator cable
- Timer cable
- Remote exposure switch cable (optional)
- Electrical cable
- RVG cable (optional)
- Connection cable to the "Ready state" lamp

External cables must pass through the openings (c) in the wall framework.




Position the insulating sheet on the 4 plastic pins (d) and on the 2 metal ground pins (e) (see next diagram).

Install the board in the wall framework on the nylon and metal pins.

Mount the nylon screws (f) on the plastic pins (d) and the metal screws (g) on the metal pins (e).


Proceed to the assembly of the main power board and its adaptor. Install on the power board the spacers and fix the adapter board for the connection of the control timer cable.






3.2.2 - Wiring

Connect all the wires according to the wiring diagrams below.



- Connect RVG link cable to the RVG connector (optional)
Yellow wire on 1
White wire on 2
Neutral blue on 3

, Connect the control timer cable to the power board connector J1 (20) on the adaptor board of the power board.



f „ Electrical wiring: MAIN PWR connector on the board and On/Off button



Main Power **f**

- Power supply cable: Ground on 1 (reference \oplus)
- Power supply cable: Neutral on 3 (reference N)
- Power supply cable: Phase on 4 (reference L)

On/Off button „

- On/off cable: wire 5 on 5
- On/off cable: wire 6 on 6
- On/off cable: wire 7 on 7
- On/off cable: wire 8 on 8


...Generator wiring: PWR ARM and CTRL ARM connectors on the board



- Arm cable: brown wire on PWR ARM-1
- Arm cable: orange wire on PWR ARM-3
- Arm cable: blue wire on PWR ARM-6

- Arm cable: purple wire on PWR ARM-7
- Arm cable: green-yellow wire on PWR ARM-8
- Arm cable: white 4-pin plug on CTRL ARM

Make one turn with the brown and orange wires around the ferrite provided with the arm.




Place the excess cable inside the extension arm.

† Remote exposure switch (optional)

Connect the 2 wires from the exposure switch link cable:

- shield in 3,
- other wire in 2.




WARNING

Follow the connection diagram. Neither of the 2 wires of the hand-switch cable should be connected to the diode terminal, represented by




‡ Warning lamp: a connection has been provided making it possible to indicate the "Ready" condition at a location away from the control position, in compliance with standards in force. To do this, connect a 110 V or 230 V bulb (depending on the operating voltage) to the terminals of connector J9 (7). Use usual means (i.e., a fuse) to protect this circuit.



4 - Installing the control timer

4.1 - Standard mounting



The control timer is ready to be plugged directly into the power board, on connector J1 (see paragraph 3.2.2).


4.2 – Option 1 - Installation with the control timer separated from the framework

In this configuration, the timer is mounted on the wall on a dedicated wall-mount assembly, at a distance from the framework up to 15 meters (16.4 yd.).



CAUTION

To be achieved, this configuration requires the addition of the remote timer box (**Cat 5154315**) to the standard configuration, that has to be ordered separately (not included in the standard kit).




Ø Description of the kit:

- 1 mounting plate with connection board
- 1 cover for the control timer
- 1 control cable measuring 15 m (16.4 yd.) to connect to the power board
- 1 link cable measuring 30 cm (11-13/16 in.) in length
- Accessories: screws, mounting ankles...

Ø Position the control cable (15 meters or 16.4 yards maximum length) between the power board and the position of the timer.

The control cable connector must be on the power board side.

- Ø Pass the connection cable through one of the dedicated holes in the base plate.
Position the mounting base on the wall and mount it through the holes made in the plate.
Use appropriate mounting media, based on the type of wall material.
Use the spirit level to check that the base is vertical.



Ø Connect each wire to the 10-pole connector one by one, according to the following color code:

- 1: orange
- 2: purple
- 3: green
- 4: black
- 5: gray
- 6: brown
- 7: yellow
- 8: red
- 9: white
- 10: blue + shield

10	9	8	7	6	5	4	3	2	1

Mount the cable properly with its clamp, and make sure that the shield is connected properly.

Ø Locate the handheld timer:

- Remove the 2 protective rubber feet (A) from the back of the remote.


- Unscrew the 2 screws (B), and remove the back cover (C).

- Remove the spiral cable (D) from the timer.

- Replace it with the supplied 30 cm (11-13/16 in.) link cable (E).

NOTE


Depending on the configuration of the exposure hand-switches, you might need to reconfigure the dipswitches on the timer. Refer to Section 5.



Ø Pass the link cable through the hole in the cover and place the front panel of the timer on the cover.

Ø Mount the front panel to the cover with the T45 screws (H) (supplied).

Ø Position the assembly (F+G) in front of the plate (I), and connect the link cable (E) to the connection board. Mount the assembly on the mounting plate with the 2 supplied screws M5-16 (K). Cover the fixing screws with the triangular covers (L).



4.3 – Option 2 - Installation with a separate exposure hand-switch

In this configuration, the exposure is activated only from the exposure switch. The exposure button on the control timer is disabled.



CAUTION

To be achieved, this configuration requires the addition of the separated exposure switch (**Cat 5154307**) to the standard configuration, that has to be ordered separately (not included in the standard kit).




Ø Description of the kit:

- 1 exposure hand-switch (3) with its connection box


Ø Position the exposure switch cable between the power board and the position of the hand-switch.

The exposure switch cable (4) should be a shield cable 2 x 0.22 mm² (2 x 24 AWG). It is not supplied in the standard nor the accessory kit.



Ø Open the connection box by unscrewing the screw (a), and separate the plate (1) from its cover (2).

-
- Ø Position the plate (1) on the wall in its final position and mount it through the mounting holes (b).
 - Ø Connect the cable to the terminals on the base plate.
 - Ø Close the connection box by screwing screw (a)
 - Ø Disable the x-ray emission from the timer by setting dipswitches SW4 and SW5 to the **ON** position.



Dipswitches SW1 to SW8 are located on the handheld timer electronic board. To change the dipswitch configuration, see Chapter 5 "Timer Configuration."

IMPORTANT

Local regulations may require disabling the x-ray emission from the timer if a separate exposure button is available. Do not forget to position dipswitches SW4 and SW5 to ON.

- Ø Connect the exposure hand-switch to J7 (10) on the power board (see Chapter 3.2.2 "Electrical installation" in this Section).


4.4 – Option 3 - Installation with 1 separate control timer and 1 separate exposure switch (mainly for U.S. market)

In this configuration, the exposure is activated by acting simultaneously on the exposure switch and on the exposure button from the control timer.



CAUTION


To be achieved, these optional configurations require the addition of the remote timer box (**Cat 5154315**) and the separated exposure switch (**Cat 5154307**) to the standard configuration. These accessories have to be ordered separately. They do not come with the standard kit.




Ø Description of the kit:

- 1 separated wall-mount control timer kit (1) (see paragraph 4.2 above)
- 1 hand-switch with its connection box (2) (see paragraph 4.3 above)

- Ø Pass the control cable (maximum 15 meters, or 16.4 yards) between the connector on the power board side and the Control Timer unit (Cable A).
- Ø Pass the exposure switch shield cable 2 x 0,22 mm² (2 x 24 AWG) between the Control Timer unit and the hand-switch location (Cable B).
- Ø Pass the 2 cables (A & B) through the dedicated holes in the mounting plate (D).
- Ø Position the mounting plate (D) on the wall and mount it through the holes provided. Use appropriate mounting media, based on the type of wall material.



- Ø Position the connection box (E) of the hand-switch on the wall and mount it as described in paragraph 4.3.
- Ø Connect the exposure switch cable (B) on the timer connection board as indicated in the diagram below.




WARNING

Follow the connection diagram. Neither of the 2 wires of the hand-switch cable should be connected to the diode terminal, represented by



- Ø Close the connection box (E) and install the hand-switch as described in paragraph 4.3.
- Ø Enable the x-ray emission simultaneously from the timer and the hand-switch by setting the dipswitch SW5 to the **ON** position.



Dipswitches SW1 to SW8 are located on the handheld timer electronic board. To change the dipswitches configuration, see Chapter 5 “Timer Configuration”.

- Ø Connect the control timer unit (F) and cover the unit (D) as indicated in paragraph 4.2.


5 – Timer Configuration



WARNING

Switch the unit OFF before any operation on the control timer remote.

- Ø Remove the 2 protective rubber covers and unscrew the 2 screws to open the control timer.
- Ø Look for the dipswitches as indicated below. The default configuration is all dipswitches set to OFF.










This factory default configuration only allows the use of the handheld control timer for x-ray emission.

- Ø Modify the configuration as appropriate, as indicated in the table below:

IMPORTANT

For Australia, pre-setting of the dipswitches is not done exiting the factory. Please modify the configuration locally during installation according to your local regulations.

Configuration of the dipswitches		Description of the configuration
ON		Factory default configuration. X-rays can only be emitted by the control timer.
ON		RVG mode is disabled.
ON		Australian mode: selection of film type 1 to 6 by 2 additional buttons.
ON		Australian mode: selection of film type 1 to 5 by 2 additional buttons.
ON		U.S. mode: X-rays can only be emitted by the control timer and a separate hand exposure switch simultaneously.
ON		X-rays can be emitted either by the control timer remote or by a separated hand exposure switch.
ON		X-rays can only be emitted by a separate hand exposure switch.

- Ø Close the control timer and replace the 2 rubber covers.

SECTION 6


SETUP

Checking the Voltage

Use the voltmeter in the 300 V AC range.

Turn the unit OFF.

Connect the voltmeter to the **L** and **N** terminals of the MAIN PWR connector (ref. 9) on the power board.



Turn the unit ON.

Set the longer time (Occlusal, type film 9, 60 kV, 7 mA).

Note the voltmeter value (no-load voltage).



Stand well behind the generator.


Make an exposure.


Note the voltmeter value during the exposure (on-load voltage).

The result of the following formula: $(V_{\text{no-load}} - V_{\text{on-load}})$ must be less than 4.5 V in 100, 110 or 130 V, and less than 9 V in 230 or 240 V (corresponding to a 3% variation).

If not, refer to “Required Electrical Specifications” in Section 4, paragraph 3.

Checking the mA

Connect a 30V range voltmeter to the test points **0 mA** (Ground ) (ref. 21 – PT25) and **RTN mA** (ref. 18 – PT17) complying with the polarities, with the plus on the **0 mA** side.



The green LEDs (ref. 3) of the power board corresponding to +12 V and +5 V must both light up.



Set the longer time (Occlusal, type film 9, 60 kV, 7 mA).
 Stand well behind the generator and make an exposure.
 Note the voltmeter value.
 This value must be between 6 V and 8V; if not, contact a technician.


The measurement is the result of the product of the current (mA) multiplied by a resistance of 1 k Ω implemented on the board.

3 - End of Installation

Place the On/Off switch in the cover recess corresponding to your orientation.

Check that all the plugs and cables are fixed securely.

Install the framework cover with 4 screws as shown below and insert the screw caps.



Clean the whole unit. Use an alcohol-based product; solvent-based products must never be used.

Follow the manufacturer's recommendations for safety precautions when using the cleaning product.

IMPORTANT

Fill out all documents concerning the installation and the warranty.
 Make sure that all manuals are left on site for the operators.

SECTION 7

PREVENTIVE MAINTENANCE

1 – Cleaning and Disinfecting

For cleaning, you are advised to use a **non-corrosive** alcohol-based product, and to avoid allowing liquids to drip into the unit.

Follow the manufacturer's recommendations for safety precautions when using the cleaning product.

The usual disinfectant products can be used, but we recommend that you protect the unit from contamination by using barriers available from dental dealers.

Follow the manufacturer's recommendations for safety precautions when using the disinfectant product.

2 - Maintenance Checks

We recommend carrying out the periodic preventive maintenance operations described below. They should be performed when the unit is first installed and then annually by a qualified technician.

IMPORTANT

If your Kodak 2200 unit does not satisfy all these checks, refer to the "Corrective Maintenance" Section in this guide, or contact your distributor. In the meantime you should not use the equipment.

2.1 - Generator

Check that the specifications label is legible and intact.

Check that there are no oil leaks.

2.2 – Mounting the complete unit to the wall

Check that all the labels are clearly legible.

Check that the wall framework is securely attached to the wall. If it is not, check the installation exactly as stipulated in the procedure described in the Section 5, "Installation."

If the installation includes a 825 mm (32-1/2 in.) extension arm, check that the 360° rotation of the scissor arm is restricted by the limitation system installed on the extension arm to prevent the power supply cable from being ripped out.

2.3 - Flexible movement

Check that the arm is flexible in all positions, and that it remains immobile when no longer handled.

2.4 - Control timer and electrical installation

Check that the symbols are always clearly legible.

Check that the handheld timer cable and the power supply cable are in good condition.

Check that the ground is correctly installed.

Check that the exposure switch returns to its initial position after use.

2.5 - Operation

Turn the timer ON.

Check that the green LED "on" is lit.



Program the exposure setting (Occlusal, type film 9, 60 kV, 7 mA).

Stand well behind the generator and take an exposure.

Check that the x-ray emission light comes on during the exposure, that the counter counts down to zero and that the audible signal stops when the x-ray emission terminates.



Make another exposure and check that when the exposure switch is released before the end of the exposure the display unit indicates "OP.ERROR" and an audible signal different from the previous one is emitted.

Self-test

- a. Turn the unit OFF.
- b. Press the RVG key while simultaneously turning the unit ON.
- c. As soon as the first light comes on, release the RVG key.

All the functions and indicators lights of the control unit will be tested one by one. The audible alarm and display unit are also tested. At the end of this test, the number of exposures carried out since the unit was put into operation is displayed.

A short beep indicates the that test has ended.

IMPORTANT

If the result of one of these checks is not satisfactory, we recommend that you contact an authorized technician for assistance. In the meantime, do not use the equipment.

At the end of this verification, make sure that the User's Guide remains with the unit.


SECTION 8

CORRECTIVE MAINTENANCE

1 - Main Components


Connectors

MAIN PWR	(9)	Electrical power supply and On/Off button
PWR ARM	(19)	Generator cable
CTRL ARM	(19)	Generator cable (safety devices)
RVG	(11)	Cable for RVG
J1	(20)	Timer
J9	(7)	Ready state lamp connector




(10)	Remote exposure switch cable
------	------------------------------

Jumpers




(1)	Adaptation to 110 V voltage
-----	-----------------------------

or



Adaptation to 230 V voltage



(2)	Strap if 110V only
-----	--------------------

LEDS

D48	(3)	+12 V
D51	(3)	+5 V
D41	(17)	Heating voltage
D65	(15)	12V non isolated

Test Points

PT10	(15)	+5 V non isolated
PT11, PT12	(16)	H Bridge control
PT13	(5)	KV return
PT14	(14)	KV ref
PT15	(5)	IHEAT return
PT16	(17)	VHEAT return
PT17	(18)	MA return
PT18	(18)	Relay R control
PT19	(3)	+12 V
PT20	(3)	+5 V
PT21	(13)	Vdac ref
PT22, PT23, PT24	(12)	I2C bus
PT25	(21)	Ground – 0 mA

Fuses


Protection of power components

F1 (8) in 230 V/250 V: 5 A 5 x 20 fast,
with high breaking capacity - UL & IEC60127-2.
Ref: 0001.1014 (SCHURTER) or 1942100002 (WICKMANN).

in 100 V/110 V/130 V: 10 A 5 x 20 fast,
with high breaking capacity - UL & IEC60127-2.
Ref: 0001.1011 (SCHURTER) or 1941500002 (WICKMANN).

Neutral line

F 2 (6) neutral links (except for mobile option).



2 - Diagnosis

2.1 - Timer self-test

Turn the unit OFF.
Press the RVG key while simultaneously turning the unit ON.
As soon as the first light comes on, release the RVG key.

All the functions and indicators lights of the control unit will be tested one by one. The audible alarm and display unit are also tested. At the end of this test, the number of exposures carried out since the unit was put into operation is displayed.

A short beep indicates that the test has ended.

If the test is faulty, the corresponding error message will be displayed.

IMPORTANT

If the result of one of these checks is not satisfactory, we recommend that you contact an authorized technician for assistance. In the meantime, do not use the equipment.

2.2 - Error messages

ERROR MESSAGE	CAUSE	HOW TO CANCEL
COOLING	Cooling cycle; this message can appear during a period of intensive use.	Do not turn the system off. The error message will disappear when the system returns to a satisfactory temperature.
CAUTION: If you turn off power to the system, the microprocessor doesn't calculate the cooling time, and for safety reasons considers that the system has not gone into the cooling cycle.		
OP. ERROR Plus audible alarm	Release of the exposure control switch before the end of the exposure. The display shows the remaining exposure time.	Press any key on the arch to stop the alarm.
KV ERROR	The generator's high voltage value is more than 10% below the required value.	Turn the system off and restart. If the problem persists, call a qualified service technician and discontinue using the equipment.
POWER ERROR	Over- or Under-voltage of main power, power supply problem for the filament, or cut filament.	Turn the system off and restart. If the problem persists, call a qualified service technician and discontinue using the equipment.

2.3 – Trouble-shooting hints


PROBLEM	CAUSE	SOLUTION
Nothing lights up	Unit is disconnected	Connect the unit
	Fuse F1 is burned out or defective	Replace the fuse
	Main circuit breaker is OFF	Turn ON the circuit breaker
Nothing lights up on the control unit	Control unit is disconnected	Connect the control unit
	Fuse F1 is burned out or defective	Replace the fuse
	The control unit is defective	Replace the control unit
No x-ray emission	The generator is cooling	Wait for the “ COOLING ” message to disappear
	The exposure switch is defective	Replace the exposure switch
Emission OK but exposure is too light, or even white	Wrong receptor type is being used	Change the receptor type (refer to User's Manual)
	Generator is incorrectly positioned	Adjust the position of the generator
	Exposure time is too short	Increase the exposure time
	Development time is too short	Increase the development time (refer to development instructions)
	Developer is too cold	Heat the developer
	Developer is too old or diluted	Change the developer
	RVG key is incorrectly selected	Verify your exposure settings (refer to the procedure in the User's Manual)
	Receptor is facing the wrong way	Reposition the receptor
	The unit was incorrectly installed	Check the installation
Emission OK, but exposure is too dark	Wrong receptor type is being used	Change the receptor type (refer to User's Manual)
	Exposure time is too long	Decrease the exposure time
	Development time is too long	Decrease the development time (refer to development instructions)
	Developer is too hot	Cool the developer
	Developer is too concentrated	Change the developer
	RVG key is incorrectly selected	Verify your exposure settings (refer to the procedure in the User's Manual)

3 – Part Replacement Procedures


3.1 - Changing the generator's power supply cable

0 - Turn the unit OFF.


1 - Remove the wall framework's cover by removing the screw covers and the 4 fixing screws.
Take the On/Off switch out of its recess.




2 - Remove the cable-clamp from the pivot bloc.




3 - Disconnect the arm's cable from the PWR ARM and CTRL ARM connectors.



4 - Remove the two plastic cable covers on the extension arm, revealing the arm cable.



5 - Take the cable out of the extension arm and let it hang from the end of the scissor arm.




6 - Remove the scissor arm covers in the following order:

(a) Covers of the generator support and the plastic cap.


(b) Cover for the central hinge, via 2 screws accessible from underneath.

(c) The 2 covers on each part of the scissor arm.




7 - Remove the 4 cable-clamps from the scissor arm, by punching in the pins with a 2 mm pin punch.

Note: the pins will fall inside the arm. Make sure you get them in order to replace the cable clamps.




8 - Take out the cable, working toward the generator. Disconnect the generator plug and remove the cable by pulling the connector. Recover the cable-clamp pins.




H214_0007AC


9 - Take the new cable. Feed it through the top opening of the generator support. Be careful which way the connectors are positioned. The cable must not be twisted during installation. You must therefore check the position of the generator connector. It is possible to place the scissor arm horizontally to facilitate the installation of the cable. Connect the cable to the generator connector.




10 - Place the scissor arm in a vertical position. Place the bottom cable-clamp, on the generator side, using 2 pins.



Extend the arm horizontally.
Place the other 3 cable-clamps with 2 pins each.
The cable must not be slack, but it must not be taut either.




Slip the cables into the extension arm through the corresponding openings. Feed the cable through the arm's support axis. Connect the cable to the board's PWR ARM and CTRL ARM connectors. Fix the cable to the arm support using the cable-clamp. Place the cable in excess in the extension arm. Test the connection. Place the arm in different positions to check any stress on the cable.



WARNING
Do not forget to reinstall properly the ferrite (brown and orange wires)


Reassemble the various covers in the order shown in the figure below.



H214_0008HC


3.2 - Generator replacement

Remove the 2 plastic covers on each side of the arm, as well as the plastic cap.



Place the scissor arm in the vertical position.

Strap the scissor arm so that it does not suddenly spring back when the generator is removed.



H214_0010GC

H214_0011GC

H214_0012GC

3.3 - Arm stability

CAUTION

After a certain amount of time, the arm may become unstable. In this case adjust the spring on the lower arm as described below (for the front arm, refer to Section 3.5 for front jack replacement).


Remove the access plate (A) with a screwdriver.

Use the 4 mm (.16 in.) diameter metal strip (B) supplied with the unit to adjust the spring (C) by turning the nut:

In the direction shown on the drawing if the arm goes down.

In the other direction if the arm goes up.

Replace the access plate (A).




3.4 - Arm flexibility adjustment

Remove the access plate (A) with a screwdriver.

Tighten the set screw to harden the movement.
Loosen the screw to smoothen the movement.

Replace the access plate (A)



3.5 - Front jack replacement

Remove the two plastic covers (E).

Place the front arm in a vertical position.

Remove the plastic cover (F) on the front arm using a 4 x 150n (.16 x 5.9 in.) screwdriver.

Remove the access plate (A) using the screwdriver, by pressing on its base, reference P.

Remove the circlips (B) and (G).

Support the arm.
Withdraw axles (H) and (C).

Extract the jack (D).

Take the new jack and insert it rod-end first in the arm opening. Place it in position and reinsert axle (H).


Take axle (C). Place the jack end between the plates. Insert the axle on this assembly.

Place the circlips (B) and (G) back in place.

Replace the access plate (A).

Place the front arm in a vertical position again. Put the cover (F) back in position.

Install the two covers (E) back on the arm's end and the plastic cap.




APPENDIX

APPENDIX 1

DIMENSIONAL DIAGRAMS

Dimensions


Space requirements - frame oriented to the right



H214_0014HC

EXTENSION	R	A	B	C
CG645	470 mm (18-1/2 in.)	651 mm (25-5/8 in.)	622 mm (24-1/2 in.)	757 mm (29-13/16 in.)
CG646	648 mm (25-1/2 in.)	829 mm (32-5/8 in.)	800 mm (31-1/2 in.)	935 mm (36-13/16 in.)
CG648	825 mm (32-1/2 in.)	1 006 mm (39-5/8 in.)	977 mm (38-1/2 in.)	1 112 mm (43-3/4 in.)


Space requirements - frame oriented to the left



H214_0015HC

EXTENSION	R	A	B	C
CG645	470 mm (18-1/2 in.)	651 mm (25-5/8 in.)	622 mm (24-1/2 in.)	757 mm (29-13/16 in.)
CG646	648 mm (25-1/2 in.)	829 mm (32-5/8 in.)	800 mm (31-1/2 in.)	935 mm (36-13/16 in.)
CG648	825 mm (32-1/2 in.)	1 006 mm (39-5/8 in.)	977 mm (38-1/2 in.)	1 112 mm (43-3/4 in.)


Space requirements - frame oriented vertically



H214_0016HC

EXTENSION	R	A	B	C
CG645	470 mm (18-1/2 in.)	666 mm (26-1/4 in.)	637 mm (25-1/16 in.)	772 mm (30-3/8 in.)
CG646	648 mm (25-1/2 in.)	844 mm (33-1/4 in.)	815 mm (32-1/16 in.)	950 mm (37-3/8 in.)
CG648	825 mm (32-1/2 in.)	1 021 mm (40-3/16 in.)	992 mm (39-1/16 in.)	1 127 mm (44-3/8 in.)


Space requirements - side view



H214_0017GC


EXTENSION	R	A	B	C
CG645	470 mm (18-1/2 in.)	615 mm (24-1/4 in.)	586 mm (23-1/16 in.)	721 mm (28-3/8 in.)
CG646	648 mm (25-1/2 in.)	793 mm (31-1/4 in.)	764 mm (30-1/16 in.)	899 mm (35-3/8 in.)
CG648	825 mm (32-1/2 in.)	970 mm (38-1/4 in.)	941 mm (37-1/16 in.)	1 076 mm (42-3/8 in.)

Clearance space - side view



EXTENSION	R	A	B
CG645	470 mm (18-1/2 in.)	1 700 mm (67 in.)	1 047 mm (41-1/4 in.)
CG646	648 mm (25-1/2 in.)	1 880 mm (74 in.)	1 225 mm (48-1/4 in.)
CG648	825 mm (32-1/2 in.)	2 050 mm (80-3/4 in.)	1 402 mm (55-1/4 in.)


Clearance space - bottom view



H214_0019HC

EXTENSION	R	A	B	C	D
CG645	470 mm (18-1/2 in.)	1 700 mm (67 in.)	1 383 mm (54-7/16 in.)	785 mm (30-15/16 in.)	3 257 mm (128-1/4 in.)
CG646	648 mm (25-1/2 in.)	1 880 mm (74 in.)	1 758 mm (69-1/4 in.)	963 mm (37-15/16 in.)	3 617 mm (142-7/16 in.)
CG648	825 mm (32-1/2 in.)	2 050 mm (80-3/4 in.)	2 126 mm (83-11/16 in.)	1 140 mm (44-7/8 in.)	3 975 mm (156-1/2 in.)


Installation for two seats



H214_0020DC

EXTENSION	R	L
CG645	470 mm (18-1/2 in.)	821 mm (32-5/16 in.)
CG646	648 mm (25-1/2 in.)	1 000 mm (39-3/8 in.)
CG648	825 mm (32-1/2 in.)	1 176 mm (46-5/16 in.)

Installation for pass through (CL045) (for U.S. only)



H214_0021DC

EXTENSION	R	L
CG645	470 mm (18-1/2 in.)	821 mm (32-5/16 in.)
CG646	648 mm (25-1/2 in.)	1 000 mm (39-3/8 in.)
CG648	825 mm (32-1/2 in.)	1 176 mm (46-5/16 in.)

APPENDIX 2

Technical Specifications

1 - Technical characteristics according to IEC standard 60601-2-7 (2002)

Manufacturer

TROPHY
4, rue F. Pelloutier - Croissy-Beaubourg
77435 Marne-la Vallée Cedex 2 (France)

For :

Carestream Health, Inc.
150 Veronal Street
Rochester
New-York – USA 14608

Models

Dental x-ray diagnosis devices, class 1, type B, intermittent use.

The Kodak 2200-TR system is equipped with tube TRX 708 from Trophy.
The Kodak 2200-C system is equipped with tube OCX / 65-G from CEI.

Electric power supply (during exposure)

230 - 240 V AC ($\pm 10\%$), 50 Hz, 5 A, apparent resistance 0,5 Ω
100 - 110 - 130 V AC ($\pm 10\%$), 50/60 Hz, 12 A, apparent resistance 0,2 Ω

Electric power supply (no exposure)

230 - 240 V AC ($\pm 10\%$), 50 Hz, 100mA
100 - 110 - 130 V AC ($\pm 10\%$), 50/60 Hz, 100mA

Rated high voltage and maximum corresponding current

Film mode	70 kV	7 mA
RVG mode	70 kV	4 mA

Current/voltage combinations for a maximum output power of

490 W in film mode	70 kV / 7 mA
280 W in RVG mode	70 kV / 4 mA

Rated power for exposure time of 0.1 s

Film mode	490 W
RVG mode	280 W

Rate of use

At 70 kV, 7 mA and 0.1 s and at the maximum tank temperature: one exposure every 8 seconds.

Minimum value of the current/time product in the range of conformity

0.07 mAs at 7 mA

0.04 mAs at 4 mA

Parameters selection

Film mode	70 kV / 7 mA	and	60 kV / 7 mA
RVG mode	70 kV / 4 mA	and	60 kV / 4 mA

Area of conformity to IEC standard 60601-2-7 (2002)

Reproducibility of the emitted radiation	conform
Linearity of the emitted radiation	conform
Precision in radiography	conform

Measurement conditions

kV: Indirect measurement using a kV-peakmeter.

mAs: direct measurement in the circuit using a mAs-meter.

Exposure time: indirect measurement on the kV signal at 75% of the peak value

Storage and transport conditions

Temperature	: -10°C to 60°C
Relative humidity	: 10% to 95%
Atmospheric pressure	: 700 to 1060 hPa

Dimensions and weight

Component	Dimensions	Weight
Control unit	16 x 9 x 4 cm (6-5/16 x 3-1/2 x 1-9/16 in.)	0.4 kg
Wall framework	51.4 x 18.9 x 10.8 cm (20-1/4 x 7-7/16 x 4-1/4 in.)	4.3 kg
X-ray emitting unit	43.8 x 22.6 x 12 cm (17-1/4 x 8-15/16 x 4-3/4 in.)	4.3 kg
Scissor arm	87.3 x 13.3 x 6.3 cm (34-3/8 x 5-1/4 x 2-1/2 in.)	9 kg

Scissor arm

The scissor arm is equipped with one gas jack specially designed for this particular application.

It has been proven to function correctly after more than 400,000 cycles.


Electro-magnetic Compatibility

The Kodak 2200 intraoral x-ray system complies with the European Directive 89/336/EEC and the IEC 60601.1.2 (2001) standard.


Classification: Group 1, Class B

2 - Main characteristics of the x-ray generator

Manufacturer and Type of X-ray Tube	Trophy Type TRX 708	CEI Type OCX / 65-G
Rated high voltage	70 kV	70 kV
Rated anodic power	490 W	490 W
Maximum heat accumulated in the anode	8 700 J	10 000 J
Rated value of focal spot (IEC 60336/1993)	0.7 mm (.027")	0.7 mm (.027")
Target materials	Tungsten	Tungsten
Target slope	19°	19°
Filtration due to fixed materials	0.6 mm (.023") eq. Al	0.6 mm (.023") eq. Al




Trophy TRX708 tube
Heating and cooling curves




CEI OCX/65-G TUBE
Heating and cooling curves

Equipped x-ray generator

IEC standard 60601-2-28 (1993)	Conform
Type of protection against electric shocks	Class 1
Degree of protection against electric shocks	Type B
Rated value of inherent filtration	1.5 mm (.059") eq. Al
Rated value of additional filtration	1.0 mm (.039") eq. Al
Rated value of total filtration	2,5 mm (.098") eq. Al
Beam-limiting cone, focal spot/skin distance	20 cm (7-7/8 in.)
Maximum accumulated heat	32 500 J
Maximum continuous thermal dissipation	7 W
Leaking radiation at maximum rate during one hour of use	< 0,25 mGy
Maximum field of symmetrical radiation	6 cm diameter (2-3/8 in.)
Position and tolerances of the focal spot on the reference axis	0 mm +/- 0.5 mm (0.020")




- 1 Reference axis
- 2 Focal spot
- 3 Target angle



Heating and cooling curves of the Kodak 2200 system's tube head

3 – Position of identification labels




- 1 **WARNING: MAGNETIC AREA**
- 2 Identification of the unit
- 3 Identification of the x-ray emitting unit

APPENDIX 3

SPARE PARTS LIST

DESCRIPTION	REFERENCE	CODE
GENERATOR		
Complete 110 V generator	1	SPK22-1
Complete 230 V generator	1	SPK22-2
K2200 generator cover set	2	SPCSK22
Plastic ring retainer	3	SPHY259
Fork cover set	4	SPCG349
ARM		
Scissor arm	5	SPCG640
Scissor arm for ceiling mount		SPCG692
Scissor arm cable	6	SPCP490
Scissor arm cable for ceiling mount		SPCP748
Arm jack	7	SPJW132
Arm cover set	8	SPCG352
Arm base cover set	9	SPCG353
Extension arm cover set	10	SPHY206
Cover for ceiling mount		SPBY160
WALL FRAMEWORK		
Power board + adaptation 110V	11	SPJ22-1
Power board + adaptation 230V	11	SPJ22-2
Adaptation board CJ718 -> CJ728 K2200		SPCJ724
Wall apply without power board	12	SPCG795
Plastic insulator	13	SPDE164
On/Off cable for wall framework	14	SPCP489
Wall apply cover	15	SPCL100
CONTROL TIMER		
K2200 handheld timer	16	SPCG792
RHF Timer rear panel		SPCL098
CABLES		
Timer link cable		SPCP735
Remote timer link cable		SPCP733
Remote timer cable 15m		SPCP514
MISC.		
Screw set		SPSSK2K



H214_p02-EG

Made by:

Trophy
4 rue F. Pelloutier – Croissy-Beaubourg
77435 Marne la Vallée Cedex 2 (France)
+33 1 64 80 85 00

For:

Carestream Health, Inc.
150 Veronal Street
Rochester
New-York – USA 14608



© Carestream Health, Inc., 2007
The Kodak trade mark and trade dress are
used under license from Eastman Kodak Company
