Bilirubin Total

DMSO. Colorimetric

Quantitative determination of bilirubin IVD

Store at 2-8°C

PRINCIPLE OF THE METHOD

Bilirubin is converted to colored azobilirubin by diazotized sulfanilic acid and measured photometrically. Of the two fractions presents in serum, bilirubin-glucuronide and free bilirubin loosely bound to albumin, only the former reacts directly in aqueous solution (bilirubin direct), while free bilirubin requires solubilization with dimethylsulphoxide (DMSO) to react (bilirubin indirect). In the determination of indirect bilirubin the direct is also determined, the results correspond to total bilirubin.

The intensity of the color formed is proportional to the bilirubin concentration in the sample^{1,2,3}.

CLINICAL SIGNIFICANCE

Bilirubin is a breakdown product of hemoglobin.

It is transported from the spleen to the liver and excreted into bile.

Hyperbilirubinemia results from the increase of bilirubin concentrations in plasma.

Causes of hyperbilirubinemia:

Total bilirubin: Increase hemolysis, genetic errors, neonatal jaundice, ineffective erythropoiesis, and drugs.

Direct bilirubin: Hepatic cholestasis, genetic errors, hepatocellular damage^{1,6,7}.

Clinical diagnosis should not be made on a single test result; it should integrate clinical and other laboratory data.

REAGENTS

R 1	Sulfanilic acid	30 mmol/L
	Hydrochloric acid (HCl)	50 mmol/L
	Dimethylsulphoxide (DMSO)	7 mol/L
R 2	Sodium nitrite	29 mmol/L

PRECAUTIONS

R1/RW: Corrosive (C):R35:Causes severe burns.

S26: In case of contact with eyes, rinse immediately with plenty of water and seek medical advice.

PREPARATION

Pipette 1.2 mL of R2 into R1 content. Mix avoiding foam forming and it will be ready to use (WR). Do not use the reagent before 30 min. after the reagent preparation.

STORAGE AND STABILITY

All the components of the kit are stable until the expiration date on the label when stored tightly closed at 2-8°C, protected from light and contaminations prevented during their use. Do not use reagents over the expiration date.

Signs of reagent deterioration:

- Presence of particles and turbidity.
- Color development in R 2.

ADDITIONAL EQUIPMENT

- MINDRAY BS-120 / BS-200 Autoanalyser.
- General laboratory equipment.

SAMPLES

Serum or plasma, free of hemolysis¹. Protect samples from direct light.

Stability: Bilirubin is stable at 2-8°C for 4 days and 2 months at – 20°C.

MINDRAY BS-120 / BS-200E APPLICATION

PARAMETERS			
Test	BILT	R1	250
Nº	**	R2	*
Full Name	BILIRUBIN	Sample Volume	25
Standard nº		R1 Blank	
Reac. type	Endpoint	Mixed Rgt Blank	
Pri. Wavelength	546 nm	Linearity range	0.1 - 15.0
Sec. Wavelength		Linearity limit	*
Direction	Increase	Substrate limit	*
Reac. time	1-33 / -1-45	Factor	*
Incub. Time		Prozone check	*
Units	mg/dL	q1	q2
Precision	0.01/0.1	q3	q4
		PC	Abs
CALIBRATION (Cal + Rg Blk)			
Rule	One-point Linear / Two-point linear		
Sensitivity	1		
Replicates	2		
Interval (days)	0		
Difference limit			
SD			
Blank Response			
Error Limit			
Correlation coefficient			

Blank parameter must be performed in order to get good results in CALIB screen from main menu. The blank calibration is stable until **20 days**. After this period the blank parameter must be performed again in order to validate the calibration.

QUALITY CONTROL

Control sera are recommended to monitor the performance of assay procedures: SPINTROL H Normal and Pathologic (Ref. 1002120 and 1002210).

If control values are found outside the defined range, check the instrument, reagents and calibrator for problems.

Each laboratory should establish its own Quality Control scheme and corrective actions if controls do not meet the acceptable tolerances.

REFERENCE VALUES¹

Bilirubin Total Up to 1.10 mg/dL $\cong$ Up to 18.81 μ mol/L

These values are for orientation purpose; each laboratory should establish its own reference range.

BIBLIOGRAPHY

- Kaplan A et al. Bilirubin. Clin Chem The C.V. Mosby Co. St Louis. Toronto. Princeton 1984; 1238-1241. 436 and 650.
- Malloy H T. et al. The determination of bilirubin with the photoelectric colorimeter. J. Biol Chem 1937; 112, 2; 481-491.
- Martinek R. Improved micro-method for determination of serum bilirubin. Clin Chim 1966; Acta 13: 61-170.
- Young DS. Effects of drugs on Clinical Lab. Tests, 4th ed AACC Press, 1995.
- Young DS. Effects of disease on Clinical Lab. Tests, 4th ed AACC 2001.
- Burtis A et al. Tietz Textbook of Clinical Chemistry, 3rd ed AACC 1999.
- Tietz N W et al. Clinical Guide to Laboratory Tests, 3rd ed AACC 1995.

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BIBLIOGRAFIA

- Kaplan A et al. Bilirubin. Clin Chem The C.V. Mosby Co. St Louis. Toronto. Princeton 1984; 1238-1241. 436 and 650.
- Malloy H T. et al. The determination of bilirubin with the photoelectric colorimeter. J. Biol Chem 1937; 112, 2; 481-491.
- Martinek R. Improved micro-method for determination of serum bilirubin. Clin Chim 1966: Acta 13: 61-170.
- Young DS. Effects of drugs on Clinical Lab. Tests, 4th ed AACC Press, 1995.
- Young DS. Effects of disease on Clinical Lab. Tests, 4th ed AACC 2001.
- Burtis A et al. Tietz Textbook of Clinical Chemistry, 3rd ed AACC 1999.
- Tietz N W et al. Clinical Guide to Laboratory Tests, 3rd ed AACC 1995.

PRESENTACION

Ref: MI1001042	Cont.	R 1: 5 x 30 mL
		R 2: 1 x 10 mL