

SIEMENS AG
Medical Solutions

ESPREE

PCB	Part1-No	Doc	Index	Rev	Type-No	Mod	Page
Q22	1/2	n.a.	EKO	01P	000	n.a.	--

Topsheet

Digicam *

- Digicam_Ethernet
- PWR_Digicam(0:3)
- GND

2203 Printer *

- Parallel_Printer(0:25)
- PWR_Printer(0:3)

Power Distribution

- Ser_Host(0:8)
- unused_1
- unused_2
- Inroom_Space_Mouse
- Foot_Switch
- Space_Mouse
- Ser_MPSU(0:8)
- Ser_Standby
- unused_3

ACC

- PHYSIO_BEEPER(0:2)
- DOOR_CONTACT(0:2)
- SE_CTRL(1:1)
- RCCS_BCCS_CTRL(1:5)
- MH_SEL(1:3)
- RFIS_DYN_CTRL(1:7)
- PHYSIO_DATA_CAN_INCRMENT(1:6)
- SSB1_2(0:1)
- PW(0:1)
- PICKUP_BC1(2)
- PICKUP_BC2(1:2)
- PICKUP_BC2(1:2)
- CV_OV(0:1)
- CV2(0:1)
- TTX2(0:1)
- MRSIG2(1-8)
- MRSIG2(9-16)
- MRSIG2(17-24)
- MRSIG2(25-32)
- SHIM(0:10)
- MSUP_CTRL(1:8)
- SLIO_IN_OUT(0:12)
- CAN_MPSU(0:3)
- PWR_PS_RFIS(0:3)
- PWR_SHIM_IRON
- PWR_PTAB(0:3)
- PWR_Reserve(0:3)

Main Box

- PA_SWITCH
- GRAD_Z(0:4)
- GRAD_Y(0:4)
- GRAD_X(0:4)
- GRAD_Z(0:4)
- GRAD_Y(0:4)
- GRAD_X(0:4)
- N_COILERR_FOC(1:1)


GPA

Video Distribution

- VIDEO(0:1)
- GND
- Ser_Water(1:2)
- Helium(1:2)

Connections:

- <W0570> DOOR_CONTACT(0:2) <<>
- PWR_COLO_HEAD(0:3) <<>
- PWR_PTAB(0:3) <<>
- <W0190> SSB1_3(0:1) <<>
- CAB_INTERCOM(0:25) <<>
- <W1425> SE_CTRL(1:1) <<>
- <W1425> RCCS_BCCS_CTRL(1:5) <<>
- <W1420> MH_SEL(1:3) <<>
- <W1421> RFIS_DYN_CTRL(1:7) <<>
- <W1424> Physio_Data(1:6) <<>
- PW(0:1) <<>
- PICKUP_BC_1 <<>
- PICKUP_BC_2 <<>
- CV_OV(0:1) <<>
- CV2(0:1) <<>
- TTX2(0:1) <<>
- MRSIG2(1-8) <<>
- MRSIG2(9-16) <<>
- MRSIG2(17-24) <<>
- MRSIG2(25-32) <<>
- SHIM(0:10) <<>
- CAN(1:2) <<>
- FOC CONVERTER <<>
- <W0415> GRAD_I_OVERTEMP <<>
- <W0650> PWR_PS_RFIS(0:3) <<>
- <W2004> MSUP_CTRL(1:8) <<>
- <W0680> <<>
- <W001> GRAD_Z(0:4) <<>
- <W002> GRAD_Y(0:4) <<>
- <W003> GRAD_X(0:4) <<>
- <W1430> N_COILERR_FOC(1:1) <<>
- VIDEO(0:1) <<>
- GND <<>
- Ser_Water(1:2) <<>
- Helium(1:2) <<>



Topsheet Examination Room Espree

