

	FILM UNIT OP200		DIGITAL UNIT OP200 D	
kV- display Ch / Sy / Er	mA- display error number	Second / dose - display	mA- display error number	Second / dose - display
		Error name FILM UNIT		Error name DIGITAL UNIT
Sy	22	Arc TUBEHEAD ARC	22	Arc TUBEHEAD ARC
Sy	23	Inu INVERTER FAILURE	23	Inu INVERTER FAILURE
Sy	24	FIL FILAMENT FAILURE	24	FIL FILAMENT FAILURE
Sy	25	AEC AEC FREQUENCY ERROR	25	AEC AEC FREQUENCY ERROR
Sy	27	Por ROTATION MOVEMENT FAILURE	27	Por ROTATION MOVEMENT FAILURE
Sy	28	PoC CASSETTE MOVEMENT FAILURE	28	CCd CAMERA ERROR
Sy	29	PoL LINEAR MOVEMENT FAILURE	29	PoL LINEAR MOVEMENT FAILURE
Sy	30	PoH CASSETTE HOLDER MOVEMENT FAILURE	30	PoC CEPH MOVEMENT FAILURE
Sy	31	PoU Z MOVEMENT FAILURE	31	PoU Z MOVEMENT FAILURE
Sy			32	PoA CEPH BEAM ALIGNMENT FAILURE
Sy	33	nSy NO SYNC ERROR	33	nSy NO SYNC ERROR
Sy			34	